

## INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME L'ABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80 - 006 498	Instrument Location OCSO	
Reason for Install/Repair		
Install After Receiving From Crime Laboratory 🔲 Install After	er Location Change	
Other (Specify)		
Check When Done:		
1. Surge Protector Installed/Property Grounded.		
2. Telephone Line Connected to Intoxilyzer® 8000.		
3. Breath Tube Heated.		
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Locati	on; Level 2, Function E).	
5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
6. Run Tests:		
A. Print Test (Level1, Function P).		
B. ACA Test (Level 1, Function C).		
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	iring Test).
7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.
9. File Previous Intoxilyzer® Record (SFN504096, Form 120-0	· · · · · · · · · · · · · · · · · · ·	
10. Send the Following to the Crime Laboratory: Completed Into 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair C	theckout (SFN59281, Form
Field Inspector Signature BUH II and		Date 10-10-28
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on ficulty of Burleigh, North Dakota, is certified to be a true and correct co	ile at the Office of Attorney General, Crir	

FORM 104-G Uploaded 27Oct2025

Anna Narehood

Reviewed/Certified By

Certified Date

27Oct2025

Intoxilyzer Test Record and Checklist. NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006498
Location = DUNN 8164.14.00 09/16
10/10/2025 10:37

\*\*\*\*\*\*\*\*\* Printer Test \*\*\*\*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&\*() +?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No · Display Volume? No. Disable On Memfull? Yes # of Print Copies? Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 01524080A4 Standard Cyl #? 21 Standard Expiration? 03/05/2026 Oper No? 131919

Flow Cal. Date: 08/18/2015 Slope 669 Intercept -637469

IR Calibration Date: 08/18/2015
3um 9um
0th Coef(\*100): -28161 -22508

1st Coef(\*100): 286657 143843 2nd Coef(\*100): 2341 1281 H2O adj(mg/l\*10k): 776 451

\*\*\*\*\*\* Printer Test End \*\*\*\*\*\*\*

Operator Signature BETH ISAAK

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006498 Location = DUNN 8164.14.00 09/16. 10/10/2025

Alcohol Analyzer

## DRY CAL CHECK

Time
10:38
10:39
10:39
10:40
10:40
10:41
10:41

Lot No = 01524080A4

Cyl No = 21

Exp Date = 03/05/2026

County = 13

Oper No. = 131919

Operator Signature BETH ISAAK

Remarks:

Form 106-I8000

. Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 SN 80-006498 Location = DUNN 8164.14.00 09/16 10/10/2025

Alcohol Analyzer 8164.14.00 09/16 12:33

Test	AC	Time
01 Room Air	0.000	12:36
02 *Subject Test	RFI*	12:36
03 Room Air	RFI*	12:37

\*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 01/01/2001

Sub Sex = Unknown Weight = NA

Test = OTH

Cit = NA

Dr. Lic. = ND/NA

Lot No = 01524080A4

Cyl No = 21

Expiration Date = 03/05/2026

County = 13

Oper No. = 131919

Operator Signature BETH ISAAK

Remarks:

Form 106-I8000