

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number	Instrument Location	
80-005937	Belcourt PD	
Reason for Install/Repair		
	er Location Change	
Other (Specify)		
Check When Done:		
☑ 3. Breath Tube Heated.		
✓ 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).		
✓ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).		
★ 6. Run Tests:		
⋈ B. ACA Test (Level 1, Function C).		
☑ C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).		
▼ 7. Repair and/or Maintenance Performed (if any):		
8. Complete the Top Portion of the Intoxilyzer® Record (SFN50)	0496, Form 120-G) and Place it by the li	ntoxilyzer® for Use.
☑ 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G)) at the Intoxilyzer® Location at the Age	ncy.
10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.		
Field Inspector Signature		Date
Mild Det		6/17/2025
Crime Laboratory Use Only		
This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from		
the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.		
Reviewed/Certified By Anna Ingemansen	NN ligeningen	Certified Date 20Nov2025
Ailia iliyellialiseli (///		ZUNOVZUZO

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005937
Location = BELC 8164.14.00 09/16
06/17/2025 20:37

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled Start Test Sequence: DABACABA Display Prelim Rslt? Yes Display Third Digit? Yes Inhib Printer (Y/N)? No Display Volume? No Disable On Memfull? Yes # of Print Copies? 1 Select Std (D/W/I)? Dry Standard Value? 0.080 Standard Lot #? 08425080A1 Standard Cyl #? 20 Standard Expiration? 05/05/2027 Oper No? 130918

Flow Cal. Date: 11/17/2017 Slope 606 Intercept -534563

IR Calibration Date: 11/17/2017
3um 9um

Coef(*100):	-30534	-23772
Coef(*100):	294336	142357
Coef(*100):	-906	850
adj(mg/l*10k):	205	250
	Coef(*100): Coef(*100): Coef(*100): adj(mg/l*10k):	Coef(*100): 294336

******* Printer Test End *******

Operator Signature MITCH SLATER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = BELC
06/17/2025

Alcohol Analyzer SN 80-005937 8164.14.00 09/16 20:37

DRY CAL CHECK

$T\epsilon$	est		AC	Time
01	Room	Air	0.000	20:38
02	Std.	Gas	0.079	20:38
03	Room	Air	0.000	20:39
04	Std.	Gas	0.078	20:39
05	Room	Air	0.000	20:40
06	Std.	Gas	0.078	20:40
07	Room	Air	0.000	20:41

Lot No = 08425080A1

Cyl No = 20

Exp Date = 05/05/2027

County = 63

Oper No. = 130918

Operator Signature MITCH SLATER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer	Alcohol Analyzer
North Dakota Model 8000	SN 80-005937
Location = BELC	8164.14.00 09/16
06/17/2025	20:41

Test	AC	Time
01 Room Air	0.000	20:43
02 *Subject Test	RFI*	20:43
03 Room Air	0.000	20:43

*Invalid Test Inhibited - RFI

Sub Name = TEST, SUBJECT A

Sub DOB = 01/10/1999

Sub Sex = Male Weight = 250

Test = DUI
Dr. Lic. = ND/TEST

Lot No = 08425080A1

Cyl No = 20

Expiration Date = 05/05/2027

County = 63

Oper No. = 130918

Cit = 061725

Operator Signature MITCH SLATER

Remarks:

Form 106-I8000