

INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59281 (06/2018)

Serial Number 80-005939	Instrument Location Williston				
Reason for Install/Repair					
☐ Install After Receiving From Crime Laboratory ☐ Install Afte	er Location Change				
X Other (Specify) Transported to Watford City for Update by Crime Laboratory					
Check When Done:					
☑ 2. Telephone Line Connected to Intoxilyzer® 8000.					
☑ 3. Breath Tube Heated.					
4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).					
☑ 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).					
☑ 6. Run Tests:					
☑ A. Print Test (Level1, Function P).					
☑ B. ACA Test (Level 1, Function C).					
C. Radio Frequency Interference (RFI) Test (CMS Mode or	Level 1, Function B or C; Key Radio Du	uring Test).			
☑ 7. Repair and/or Maintenance Performed (if any): ✓ A					
8. Complete the Top Portion of the Intoxilyzer® Record (SFN5	0496, Form 120-G) and Place it by the I	ntoxilyzer® for Use.			
	6) at the Intoxilyzer® Location at the Age	ency.			
10. Send the Following to the Crime Laboratory: Completed Int 104-G), Print Test, ACA Test, and RFI Test.	oxilyzer® 8000 Installation and Repair 0	Checkout (SFN59281, Form			
Field Inspector Signature		Date			
15 2		11/10/2025			
Crime Laboratory Use Only					
This installation has been reviewed and the instrument is approved to be the date the Field Inspector performed the installation. This record on for County of Burleigh, North Dakota, is certified to be a true and correct county.	ile at the Office of Attorney General, Cri				
Reviewed/Certified By Anna Ingemansen	N ligeningen	Certified Date 12Nov2025			

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005939
Location = WILL 8164.14.00 09/16
11/10/2025 14:11

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-= |
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

abcdefghijklmnopqrstuvwxyz1234567890-=| ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*() +?

Current Instrument Setup

Data Entry Mode:

Start Test Sequence:
Display Prelim Rslt?
Display Third Digit?
Inhib Printer(Y/N)?
Display Volume?
No
Disable On Memfull?
of Print Copies?
Select Std (D/W/I)?
Standard Value?
Standard Lot #?
Standard Cyl #?
Standard Expiration?
Oper No?

Enabled
DABACABA
Yes
DABACABA
Yes
Inhib Printer(Y/N)?
No
Display Volume?
No
No
Display Volume?
No
Oper No
Op

Flow Cal. Date: 06/18/2020 Slope 611 Intercept -592495

IR Calibration Date: 06/18/2020

3um 9um

Oth Coef(*100): -15726 -19510

1st Coef(*100): 296627 143165

2nd Coef(*100): -369 494

H2O adj(mg/l*10k): 220 240

****** Printer Test End *******

Operator Signature PRESTON LANGER

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WILL
11/10/2025

Alcohol Analyzer SN 80-005939 8164.14.00 09/16 14:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:14
02 Std. Gas	0.080	14:15
03 Room Air	0.000	14:15
04 Std. Gas	0.079	14:16
05 Room Air	0.000	14:16
06 Std. Gas	0.081	14:16
07 Room Air	0.000	14:17

Lot No = 01524080A4

Cyl No = 37

Exp Date = 03/05/2026

County = 53

Oper No. = 130996

Operator Signature PRESTON LANGER

276

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-005939
Location = WILL 8164.14.00 09/16
11/10/2025 14:18

Test			AC	Time
01	Room	Air	RFI*	14:19
02	Room	Air	RFI*	14:19

*Invalid Test Inhibited - RFI

Sub Name = RFI TEST, RFI TEST RFI TEST

Sub DOB = 01/01/1900

Sub Sex = Male

Weight = NA

Test = DUI

Cit = NA

Dr. Lic. = ND/NA

Lot No = 01524080A4

Cyl No = 37

Expiration Date = 03/05/2026

County = 53

Oper No. = 130996

Operator Signature PRESTON LANGER

Remarks:

Form 106-I8000