

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION

INTOXILYZER® 8000 CALIBRATION

Intoxil	zer® 8000 Serial Number: 80-00 6000 Calibration Location: TOXL
A.	Pre-Calibration Download Data Upload Operator File Buttery Check
	Was the external battery pack replaced? Yes or No
	Replaced Simulator O-Ring Replaced Breath Tube O-Ring Yes or No Yes or No
B.	Diagnostics passed and instrument in "Ready" mode Breath tube heated Date, time and location code (Level 2,E). Re-set if necessary. Time Zone: CST or MDT (Time on test records will be in time zone circled). Print test (Level 1,P). Sign and attach test record. Tank monitor (Level 3,D,G). Display: psi Regulator: psi Display and Regulator ± 50 psi of each other Yes or No Gas tank tare necessary? Yes of No If Yes, display readings after tare (Level 3,M,C,G): Display: psi Regulator: psi
C.	Tests (Sign and attach test records): Configure simulator for the following test (Level 1,S). Wet Calibration Check - Low AC (Level 1,C) Known Value ≤ 0.03 AC: 0.020 AC Sim. Ser #: MP 5289 Lot #: 2023110 Exp. Date: 112825 Results ± 0.005 of known AC

Intoxilyzer 8000 Calibration

Document ID: 11871 Revision: 3

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Status: Published

Approved By: Laboratory Director

Date Approved: 02/29/2024 Page 1 of 2

UNCONTROLLED WHEN PRINTED

2.	Configure simulator for the following test (Level 1,S). Wet Calibration Check - High AC (Level 1,C) Known Value ≥ 0.25 AC: 0.400 AC Sim. Ser #: MP\u037 Lot #: 2023\20 Exp. Date: 12\11\25 Results ± 5% AC of known AC				
3.	Configure dry gas standard for the remaining tests (Level 1,S). Known Value: 0.080 AC Gas Cylinder Lot #: 198 2508 0A2 Cylinder #: 20 Exp. Date: 915127				
4.	Enterferent Check (Level 1,B) Known Value: 0.10 AC + 0.05% Acetone Sim. Ser #: DRTIII Lot #: 1CS 9 Exp. Date: NIA Display reads "Interferent Detect" RFI Check (CMS Mode)				
5.	☐ Display reads "RFI Detect"				
6.	Test 1 0.080 Test 4 0.080 Test 7 0.080 Test 8 0.080 Test 8 0.080 Test 8 0.080 Test 9 0.080 Test 9 0.080 Average 0.080 Average 0.080				
D. Rem	arks/Maintenance: NIA				
Instrument is acceptable to be used in the field. Yes or No If No, state reason(s) why:					
Calibration Analyst Signature Date Anull Prikabilin Reviewer Date					

Intoxilyzer 8000 Calibration

Laboratory Unit: Toxicology Unit - Breath Alcohol Section

Approved By: Laboratory Director UNCONTROLLED WHEN PRINTED

Document ID: 11871 Revision: 3

Status: Published

Date Approved: 02/29/2024 Page 2 of 2

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
11/21/2025 08:43

******** Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup

Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes Display Third Digit? No Inhib Printer(Y/N)? Display Volume? No Yes Disable On Memfull? # of Print Copies? 1 Select Std (D/W/I)? Dry 0.080 Standard Value? Standard Lot #? 19825080A2 Standard Cyl #? 20 09/05/2027 Standard Expiration? 133237 Oper No?

Flow Cal. Date: 05/16/2023 Slope 686 Intercept -761238

IR Calibration Date: 11/20/2025
3um 9un

Oth Coef(*100): -17770 -27082 1st Coef(*100): 265146 136552 2nd Coef(*100): 2379 1305 H2O adj(mg/l*10k): 412 507

******* Printer Test End ********

Operator Signature

Remarks: Print Test

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = TOXL 11/21/2025

Alcohol Analyzer SN 80-006666 8164.14.00 09/16 08:46

WET CAL CHECK

Test		AC		Time
01	Room	Air	0.000	08:46
02	Std.	Sol.	0.021	08:47
	Room		0.000	08:48
04	Std.	Sol.	0.020	08:49
05	Room	Air	0.000	08:49
06	Std.	Sol.	0.020	08:50
07	Room	Air	0.000	08:50

08 Sim Temp = 34.0°C

Simul Ser No = MP5289 Std Sol No = 202311D

County = 08

Oper No. = 133237

Operator Signature ANNA INGEMANSEN

LOW AC Check-0.020AC Remarks:

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
11/21/2025

SN 80-006666 8164.14.00 09/16 09:00

WET CAL CHECK

Test		AC	
01	Room Air	0.000	09:00
	Std. Sol.	0.399	09:01
03	Room Air	0.000	09:02
	Std. Sol.	0.401	09:03
	Room Air	0.000	09:03
06	Std. Sol.	0.400	09:04
07	Room Air	0.000	09:04

08 Sim Temp = 34.0°C

Simul Ser No = MP6037Std Sol No = 202312D

County = 08

Oper No. = 133237

Operator Signature ANNA INGEMANSEN

Remarks:

High AC Check - 0.400AC

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
11/21/2025 09:42

Test	AC	Time 09:43
01 Room Air	0.000	
02 *Subject Test	INT*	09:43
03 Room Air	0.000	09:44

^{*}Invalid Test

Interferent Detected

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA Test = OTH Cit = INTERFERENT CK

Dr. Lic. = ND/TES989643

Lot No = 19825080A2

Cyl No = 20

Expiration Date = 09/05/2027

County = 08 Oper No. = 133237

Operator Signature ANNA INGEMANSEN

Remarks: Interferent check 0.10AC+0.05°10 Acetone

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006666
Location = TOXL 8164.14.00 09/16
11/21/2025 09:45

Test	AC	Time
01 Diagnostic	OK	09:45
02 Room Air	0.000	09:46
03 *Subject Test	RFI*	09:46
04 Room Air	0.000	09:47

*Invalid Test Inhibited - RFI

Sub Name = TEST, DONOR2 NONE

Sub DOB = 07/25/1998

Sub Sex = Male Weight = NA Test = OTH Cit = RFI CK

Dr. Lic. = ND/TES989643

Lot No = 19825080A2

Cyl No = 20 Expiration Date = 09/05/2027

County = 08 Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Operator Signature ANNA INGEMANSEN

Remarks:

RFI Check

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
11/21/2025

Alcohol Analyzer SN 80-006666 8164.14.00 09/16 09:48

DRY CAL CHECK

Test		AC		Time
01	Room	Air	0.000	09:49
	Std.		0.080	09:49
	Room		0.000	09:50
	Std.		0.080	09:50
05	Room	Air	0.000	09:51
06	Std.	Gas	0.080	09:51
07	Room	Air	0.000	09:52

Lot No = 19825080A2

Cyl No = 20

Exp Date = 09/05/2027

County = 08

Oper No. = 133237

Operator Signature ANNA INGEMANSEN

Remarks:

Dry Cal check #1-3

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
11/21/2025

SN 80-006666 8164.14.00 09/16 09:54

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	09:54
02	Std.	Gas	0.080	09:54
03	Room	Air	0.000	09:55
04	Std.	Gas	0.080	09:55
05	Room	Air	0.000	09:56
06	Std.	Gas	0.080	09:56
07	Room	Air	0.000	09:57

Lot No = 19825080A2

Cyl No = 20

Exp Date = 09/05/2027

County = 08

Oper No. = 133237

Operator Signature ANNA INGEMANSEN

Remarks:

Drycal check # 4-4

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = TOXL
11/21/2025

SN 80-006666 8164.14.00 09/16 09:58

DRY CAL CHECK

Test			AC	Time
01	Room	Air	0.000	09:59
02	Std.	Gas	0.080	09:59
03	Room	Air	0.000	09:59
	Std.		0.080	10:00
05	Room	Air	0.000	10:00
	Std.		0.080	10:01
07	Room	Air	0.000	10:01

Lot No = 19825080A2

Cyl No = 20

Exp Date = 09/05/2027

County = 08

Oper No. = 133237

Operator Signature ANNA INGEMANSEN

Remarks: