



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number	80-006512	Instrument Location	GFSO (SFCOC)
Reason for Install/Repair			
<input type="checkbox"/> Install After Receiving From Crime Laboratory <input checked="" type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u>Took to UND for Download, then reinstalled back at SFCOC</u>			

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): Download and New standard installed
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature	<i>[Signature]</i> #642	Date	03-05-2026
---------------------------	-------------------------	------	------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By	Anna Ingemansen	<i>[Signature]</i>	Certified Date	09Mar2026
-----------------------	-----------------	--------------------	----------------	-----------

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer           Alcohol Analyzer  
North Dakota Model 8000       SN 80-006512  
Location = GFSO                8164.14.00 09/16  
03/05/2026                      15:13

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?         No  
Display Volume?             No  
Disable On Memfull?         Yes  
# of Print Copies?          1  
Select Std (D/W/I)?         Dry  
Standard Value?             0.080  
Standard Lot #?             264437  
Standard Cyl #?             120  
Standard Expiration?        01/29/2028  
Oper No?                    133510

Flow Cal. Date:             08/20/2015  
Slope                       680  
Intercept                  -660856

IR Calibration Date:        08/21/2015  
                              3um                   9um

-----  
0th Coef(\*100):       -23051               -19261  
1st Coef(\*100):       266836               137135  
2nd Coef(\*100):       3974                1466  
H2O adj(mg/l\*10k):   785                487

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
ADAM COLE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006512  
Location = GFSO      8164.14.00 09/16  
03/05/2026      15:14

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:15
02 Std. Gas	0.081	15:15
03 Room Air	0.000	15:16
04 Std. Gas	0.082	15:16
05 Room Air	0.000	15:17
06 Std. Gas	0.081	15:17
07 Room Air	0.000	15:18

Lot No = 264437  
Cyl No = 120  
Exp Date = 01/29/2028  
County = 18      Oper No. = 133510



Operator Signature  
ADAM COLE

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006512  
Location = GFSO      8164.14.00 09/16  
03/05/2026      15:18

---

Test	AC	Time
01 Room Air	RFI*	15:20
02 Room Air	RFI*	15:20

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, TEST NA  
Sub DOB = 01/01/2000  
Sub Sex = Male      Weight = 150  
Test = OTH      Cit = 5  
Dr. Lic. = ND/NA  
Lot No = 264437  
Cyl No = 120  
Expiration Date = 01/29/2028  
County = 18      Oper No. = 133510



Operator Signature  
ADAM COLE

Remarks:

Form 106-I8000