

**NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION**

**INTOXILYZER® 8000 CALIBRATION ADJUSTMENT**

Intoxilyzer® 8000 Serial Number: 80-00 7097 Calibration Adjustment Location: TOXL

A. Pre-Adjustment

Replaced Simulator Return O-Ring Yes or No

B. Calibration Adjustment (Level 3,M,C,O)

1.  Autocalibration Printout Attached
  - Max Power Res Value  $\geq 10$
  - Auto Range Res Value  $\geq 4$
2. Simulator Solutions for Calibration Adjustment

Soln.	g/210 L	Lot No.	Exp. Date	Simulator SN
1	0.000	NA-Milli-Q H <sub>2</sub> O	NA-Milli-Q H <sub>2</sub> O	MP3066
2	0.040	<del>202401</del> <sup>2024100</sup> <sub>aci 3/21/26</sub>	22Oct26	MP3071
3	0.080	202501A	15Jan27	MP3061
4	0.100	202408F	28Aug26	MP3062
5	0.300	202408H	29Aug28	MP3058

3. 0.080 AC Calibration Gas for H<sub>2</sub>O Adjustment  
Lot No. 19825080A2 Cyl No. 20 Exp. Date: 9/5/27


4. Atmospheric Pressure
  - Displayed by Intoxilyzer® 8000 957 mbar
  - Adjusted to using barometer 958 mbar
  - Auto Calibration Report printout 957 mbar
  - Barometer Model 10510-922
  - Barometer Serial Number 250063141
  - Barometer Calibration Expiration Date 04Feb27

5.  Screen displayed "Calibration Success"
6.  Calibration Adjustment Printout Attached
  - Solution 1 Avg % Abs  $\leq 0.2500$
  - Solution 2-5 REL STD DEV  $\leq 3.000$
  - Residual (g/210 L) values for solutions 1 - 5  $\leq 0.0020$  for 3  $\mu$ m and 9  $\mu$ m channels


Dry Gas H<sub>2</sub>O adjustment sum for 3 μm and 9 μm channels within ± 10  
3 μm 3182 (Ave.) + 627 (H<sub>2</sub>O Adj.) = 3809  
9 μm 3256 (Ave.) + 553 (H<sub>2</sub>O Adj.) = 3812

C. Is an Annual Inspection due for this instrument? Yes or **(No)**  
If Yes, complete Intoxilyzer 8000 Annual Inspection (Document ID: 11698)  
If No, complete Intoxilyzer 8000 Calibration (Document ID: 11871).

Remarks/Notes: NIA  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

  
Breath Alcohol Analyst Signature

25 March 2026  
Date

  
Reviewer Signature

31 Mar 2026  
Date

TOXL  
Intoxilyzer - Alcohol Analyzer  
Model 8000 SN 80-007097  
03/24/2026 13:48:57

Auto Calibration  
Max Power Res Value = 58  
Auto Range Res Value = 40

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007097  
 03/24/2026 13:48:57

Auto Calibration

pg 1 of 2

```

    <<<<<      3um      >>>>>          <<<<<      9um      >>>>>
-----
Solution = 0.000 g/210L or 0.0000 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1  0.0900      (-0.0080)      0.1920      (0.0110)
Sample #2  0.1120      (0.0160)       0.1960      (0.0220)
Sample #3  0.1300      (0.0140)       0.2090      (0.0230)
Sample #4  0.1070      (0.0260)       0.2000      (0.0160)
Avg % Abs  0.1163      (0.0187)       0.2017      (0.0203)
STD DEV    0.0121      (0.0064)       0.0067      (0.0038)
REL STD DEV 10.398      (34.442)       3.302       (18.619)
  
```

```

-----
Solution = 0.040 g/210L or 0.1905 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1  0.8290      (0.0070)       1.5060      (0.0000)
Sample #2  0.8710      (0.0210)       1.5540      (-0.0250)
Sample #3  0.8750      (0.0330)       1.5510      (-0.0140)
Sample #4  0.8720      (0.0440)       1.5630      (0.0000)
Avg % Abs  0.8727      (0.0327)       1.5560      (-0.0130)
STD DEV    0.0021      (0.0115)       0.0062      (0.0125)
REL STD DEV 0.239       (35.215)       0.401       (96.384)
  
```

```

-----
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1  1.6320      (0.0010)       2.8140      (0.0190)
Sample #2  1.6440      (0.0130)       2.8770      (-0.0080)
Sample #3  1.5920      (0.0570)       2.8370      (0.0190)
Sample #4  1.5720      (0.0770)       2.8430      (0.0180)
Avg % Abs  1.6027      (0.0490)       2.8523      (0.0097)
STD DEV    0.0372      (0.0327)       0.0216      (0.0153)
REL STD DEV 2.319       (66.819)       0.756       (158.358)
  
```

```

-----
Solution = 0.100 g/210L or 0.4762 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1  1.9430      (0.0130)       3.4690      (0.0140)
Sample #2  1.9530      (0.0400)       3.4800      (0.0150)
Sample #3  1.9840      (0.0530)       3.4820      (0.0340)
Sample #4  1.9960      (0.0540)       3.5250      (0.0250)
Avg % Abs  1.9777      (0.0490)       3.4957      (0.0247)
STD DEV    0.0222      (0.0078)       0.0254      (0.0095)
REL STD DEV 1.122       (15.939)       0.727       (38.531)
  
```

```

-----
Solution = 0.300 g/210L or 1.4286 mg/l, Samples = 4, Discarded = 1
  Sample   % Abs      (% Abs Ref)      % Abs      (% Abs Ref)
Sample #1  5.5080      (-0.0360)      9.5090      (-0.0020)
Sample #2  5.4800      (0.0000)       9.5360      (0.0100)
Sample #3  5.4650      (0.0310)       9.5210      (0.0410)
Sample #4  5.5230      (0.0140)       9.5740      (0.0110)
Avg % Abs  5.4893      (0.0150)       9.5437      (0.0207)
STD DEV    0.0301      (0.0155)       0.0273      (0.0176)
REL STD DEV 0.548       (103.494)      0.286       (85.240)
  
```

TOXL  
 Intoxilyzer - Alcohol Analyzer  
 Model 8000 SN 80-007097  
 03/24/2026 13:48:57

Auto Calibration

pg 2 of 2

<<<< 3um >>>>			<<<< 9um >>>>		
-----			-----		
Zero Order Coef	-294.72			-287.63	
First Order Coef	2507.83			1395.55	
Second Order Coef	27.01			13.78	
-----			-----		
Act	Fit	Residual	Act	Fit	Residual
(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)	(g/210L)
0.000	-0.000	0.0001	0.000	-0.000	0.0001
0.040	0.040	-0.0002	0.040	0.040	-0.0003
0.080	0.080	0.0003	0.080	0.080	0.0001
0.100	0.100	-0.0002	0.100	0.100	0.0001
0.300	0.300	0.0000	0.300	0.300	-0.0000
-----			-----		

<<<< 3um >>>>		<<<< 9um >>>>	
-----		-----	
Solution = 0.080 g/210L or 0.3810 mg/l, Samples = 4, Discarded = 1			
Sample			
Sample #1	3206.00	3260.00	
Sample #2	3197.00	3263.00	
Sample #3	3206.00	3256.00	
Sample #4	3145.00	3250.00	
Avg	3182.6667	3256.3333	
STD DEV	32.9292	6.5064	
REL STD DEV	1.035	0.200	
H2O adjust (mg/l*10k)	627	553	

Atmospheric Pressure = 957

\*\*\*\*\*CALIBRATION SUCCESSFUL\*\*\*\*\*