



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59281 (06/2018)

Serial Number <i>80-004189</i>	Instrument Location <i>NTPD</i>
Reason for Install/Repair <input type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change <input checked="" type="checkbox"/> Other (Specify) <u><i>Install after annual inspection</i></u>	

Check When Done:

- 1. Surge Protector Installed/Property Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000.
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. Radio Frequency Interference (RFI) Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN504096, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer® 8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

Field Inspector Signature 	Date <i>4-11-26</i>
-------------------------------	------------------------

Crime Laboratory Use Only

This installation has been reviewed and the instrument is approved to be used for the analysis of breath to determine alcohol concentration from the date the Field Inspector performed the installation. This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Reviewed/Certified By <i>Anna Ingemansen</i>	Certified Date <i>01May2026</i>
---	------------------------------------

The initial print test, ACA Test and RFI Test was the incorrect Standard. This was corrected and a new print test is attached to reflect the update. This does not affect the installation of this instrument. AEI 01May2026

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004189  
Location = NTPD                    8164.14.00 09/16  
04/11/2026                            06:41

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@# \$%^ &\* ( ) \_ + ?

abcdefghijklmnopqrstuvwxyz1234567890- = |  
ABCDEFGHIJKLMN OPQRSTUVWXYZ!@# \$%^ &\* ( ) \_ + ?

Current Instrument Setup

Data Entry Mode:                    Enabled  
Start Test Sequence:                DABACABA  
Display Prelim Rslt?                Yes  
Display Third Digit?                Yes  
Inhib Printer (Y/N)?                No  
Display Volume?                    No  
Disable On Memfull?                Yes  
# of Print Copies?                  1  
Select Std (D/W/I)?                Dry  
Standard Value?                    0.080  
Standard Lot #?                    19825080A2  
Standard Cyl #?                    20  
Standard Expiration?                09/05/2027  
Oper No?                            132358

Flow Cal. Date:                    02/02/2011  
Slope                                628  
Intercept                            -452608

IR Calibration Date:                05/19/2016

	3um	9um
0th Coef (*100):	-17767	-18761
1st Coef (*100):	263439	133606
2nd Coef (*100):	1960	1140
H2O adj (mg/l*10k):	341	362

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DEMETRIUS JAMERSON

Remarks:

*Print test time is 10 mins off*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004189  
Location = NTPD                8164.14.00 09/16  
04/11/2026                        06:47

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	06:48
02 Std. Gas	RFI*	06:48
03 Room Air	0.000	06:48

\*Invalid Test  
Inhibited - RFI

Lot No = 19825080A2  
Cyl No = 20  
Exp Date = 09/05/2027  
County = 31                      Oper No. = 132358



Operator Signature  
DEMETRIUS JAMERSON

Remarks:

*RFI test time is 10 mins off.*

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = NTPD      8164.14.00 09/16  
04/11/2026      06:41

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	06:42
02 Std. Gas	0.081	06:42
03 Room Air	0.000	06:43
04 Std. Gas	0.082	06:43
05 Room Air	0.000	06:44
06 Std. Gas	0.082	06:44
07 Room Air	0.000	06:45

Lot No = 19825080A2  
Cyl No = 20  
Exp Date = 09/05/2027  
County = 31      Oper No. = 132358

  
Operator Signature  
DEMETRIUS JAMERSON

Remarks:

~~01~~ ACA Test time is 10 mins off  
Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = NTPD      8164.14.00 09/16  
04/29/2026      07:23

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPNOPQRSTUVWXYZ!@#%&^\*()\_+?

Current Instrument Setup

Data Entry Mode:      Enabled  
Start Test Sequence:      DABACABA  
Display Prelim Rslt?      Yes  
Display Third Digit?      Yes  
Inhib Printer(Y/N)?      No  
Display Volume?      No  
Disable On Memfull?      Yes  
# of Print Copies?      1  
Select Std (D/W/I)?      Dry  
Standard Value?      0.080  
Standard Lot #?      25024080A3  
Standard Cyl #?      003  
Standard Expiration?      10/05/2026  
Oper No?      132358

Flow Cal. Date:      02/02/2011  
Slope      628  
Intercept      -452608

IR Calibration Date:      05/19/2016

	3um	9um
0th Coef(*100):	-17767	-18761
1st Coef(*100):	263439	133606
2nd Coef(*100):	1960	1140
H2O adj(mg/l*10k):	341	362

\*\*\*\*\* Printer Test End \*\*\*\*\*



Operator Signature  
DEMETRIUS JAMERSON

Remarks:

*Print test*

Form 106-I8000