



NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION

INTOXILYZER® 8000 ANNUAL INSPECTION

Intoxilyzer® 8000 Serial Number: 80-00 4206 Inspection Location: TOXL

A. Pre-Inspection

1. Items with Instrument:
  - Gas Cylinder Yes or No (If Yes, Lot # \_\_\_\_\_ Cyl. # \_\_\_\_\_)
  - Keys Yes or No
  - Power Cord Yes or No
2.  Download Data
3.  Upload Operator File
4.  Current Location Code: TOXL
5.  Battery Check
  - Was the external battery pack replaced? Yes or No
6.  O-Rings
  - Replaced Simulator O-Ring Yes or No
  - Replaced Breath Tube O-Ring Yes or No

B. General Setup and Checks:

1.  Diagnostics passed and instrument in "Ready" mode
2.  Breath tube heated
3.  Date, time and location code (Level 2,E). Re-set if necessary.
  - Time Zone: CST or MDT (Time on test records will be in time zone circled).
4.  Print test (Level 1,P). Sign and attach test record.
5.  Tank monitor (Level 3,D,G).
  - Display: 485 psi Regulator: 500 psi
  - Display and Regulator ± 50 psi of each other Yes or No
  - Gas tank tare necessary? Yes or No
  - If Yes, display readings after tare (Level 3,M,C,G):
  - Display: \_\_\_\_\_ psi Regulator: \_\_\_\_\_ psi

C. Tests (Sign and attach test records):

1.  Configure simulator for the following test (Level 1,S).
  - Wet Calibration Check - Low AC (Level 1,C)
  - Known Value ≤ 0.03 AC: 0.020 AC
  - Sim. Ser #: MP5290
  - Lot #: 202501E
  - Exp. Date: 1/28/27
  - Results ± 0.005 of known AC
2.  Configure simulator for the following test (Level 1,S).

Wet Calibration Check - High AC (Level 1,C)

Known Value  $\geq$  0.25 AC: 0.300 AC

Sim. Ser #: MP3003

Lot #: 202408H

Exp. Date: 8/29/24

Results  $\pm$  5% AC of known AC

3.  Configure dry gas standard for the remaining tests (Level 1,S).

Known Value : 0.080 AC

Gas Cylinder Lot #: 20-4437

Cylinder #: 114

Exp. Date: 1/29/28

4.  Interferent Check (Level 1,B)

Known Value: 0.10 AC + 0.05% Acetone

Sim. Ser #: DR7352

Lot #: 1059

Exp. Date: None

Display reads "Interferent Detect"

5.  Invalid Sample (Level 1,B)

Display reads "Invalid Sample"

6.  RFI Check (CMS Mode)

Display reads "RFI Detect"

7.  Dry Calibration Check (Level 1,C)

Test 1 0.079 Test 4 0.080 Test 7 0.079

Test 2 0.079 Test 5 0.080 Test 8 0.079

Test 3 0.079 Test 6 0.080 Test 9 0.079

Average 0.079

Results  $\pm$  0.005 AC of known AC

D. Remarks/Maintenance: AAA del 6/10/26 Annual inspection  
computed post calibration adjustment.

Instrument is acceptable to be used in the field.  Yes or No

If No, state reason(s) why: \_\_\_\_\_

If Yes, change location code back to A.4 unless A.4 is TOXL.

Alan Lyman  
Inspector Signature

10 June 2026  
Date

Janelle Pitschler  
Reviewer

10 June 2026  
Date

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/10/2026                            14:27

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%^&\*()\_+?

Current Instrument Setup


Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Wet  
Standard Value?            0.025  
Standard Solution #?        202509A  
Sim Serial #?              MP5317  
Oper No?                    133237

Flow Cal. Date:              06/10/2009  
Slope                        690  
Intercept                    -645883

IR Calibration Date:        06/10/2026  
                                    3um                    9um

-----  
0th Coef(\*100):            -16325                    -23215  
1st Coef(\*100):            262352                    129421  
2nd Coef(\*100):            2543                      1351  
H2O adj(mg/l\*10k):        465                        605

\*\*\*\*\* Printer Test End \*\*\*\*\*

  
\_\_\_\_\_  
Operator Signature  
ANNA INGEMANSEN

Remarks:    Print Test

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
06/10/2026      14:27

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:28
02 Std. Sol.	0.019	14:32
03 Room Air	0.000	14:32
04 Std. Sol.	0.019	14:33
05 Room Air	0.000	14:34
06 Std. Sol.	0.019	14:34
07 Room Air	0.000	14:35

08 Sim Temp = 34.0°C

Simul Ser No = MP5290  
Std Sol No = 202501E  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA INGEMANSEN

Remarks:      Low AC Check  
                         0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
06/10/2026      14:47

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:47
02 Std. Sol.	0.295	14:48
03 Room Air	0.000	14:48
04 Std. Sol.	0.296	14:49
05 Room Air	0.000	14:50
06 Std. Sol.	0.298	14:50
07 Room Air	0.000	14:51

08 Sim Temp = 34.0°C

Simul Ser No = MP3003  
Std Sol No = 202408H  
County = 08      Oper No. = 133237

  
Operator Signature  
ANNA INGEMANSEN

Remarks: High AC Check  
0.300AC

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer                      Alcohol Analyzer  
North Dakota Model 8000                      SN 80-004206  
Location = TOXL                              8164.14.00 09/16  
06/10/2026    15:08

Test	AC	Time
01 Room Air	0.000	15:09
02 *Subject Test	INT*	15:10
03 Room Air	0.000	15:10

\*Invalid Test  
Interferent Detected

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                                      Weight = NA  
Test = OTH                                      Cit = INTERFERENT  
Dr. Lic. = ND/TES989643  
Lot No = 264437  
Cyl No = 114  
Expiration Date = 01/29/2028  
County = 08                                      Oper No. = 133237



Operator Signature  
ANNA INGEMANSEN

Remarks: Interferent check  
0.10AC + 0.05% Acetone

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/10/2026                            15:11

Test	AC	Time
01 Room Air	0.000	15:11
02 Invalid Sample	X.XXX	15:11
03 Room Air	0.000	15:12

\*Invalid Test - Mouth Alcohol

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                          Cit = INVALID  
Dr. Lic. = ND/TES989643  
Lot No = 264437  
Cyl No = 114  
Expiration Date = 01/29/2028  
County = 08                          Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA INGEMANSEN

Remarks: *Invalid Sample*  
*Mouth Alcohol*

Form 106-I8000

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501


CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                8164.14.00 09/16  
06/10/2026                      15:14

Test	AC	Time
01 Diagnostic	OK	15:15
02 Room Air	0.000	15:15
03 *Subject Test	RFI*	15:16
04 Room Air	0.000	15:16

\*Invalid Test  
Inhibited - RFI

Sub Name = TEST, DONOR2 NONE  
Sub DOB = 07/25/1998  
Sub Sex = Male                      Weight = NA  
Test = OTH                            Cit = RFI  
Dr. Lic. = ND/TES989643  
Lot No = 264437  
Cyl No = 114  
Expiration Date = 01/29/2028  
County = 08                          Oper No. = 133237

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature  
ANNA INGEMANSEN

Remarks:

RFI Check

Form 106-I8000


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004206  
Location = TOXL                    8164.14.00 09/16  
06/10/2026                            15:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:21
02 Std. Gas	0.079	15:21
03 Room Air	0.000	15:22
04 Std. Gas	0.079	15:22
05 Room Air	0.000	15:23
06 Std. Gas	0.079	15:23
07 Room Air	0.000	15:24

Lot No = 264437  
Cyl No = 114  
Exp Date = 01/29/2028  
County = 08                            Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA INGEMANSEN

Remarks: *Dry cal check  
#1-3*

Form 106-I8000




Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004206  
Location = TOXL      8164.14.00 09/16  
06/10/2026      15:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:28
02 Std. Gas	0.079	15:29
03 Room Air	0.000	15:29
04 Std. Gas	0.079	15:29
05 Room Air	0.000	15:30
06 Std. Gas	0.079	15:30
07 Room Air	0.000	15:31

Lot No = 264437  
Cyl No = 114  
Exp Date = 01/29/2028  
County = 08      Oper No. = 133237

  
\_\_\_\_\_  
Operator Signature  
ANNA INGEMANSEN

Remarks: Dry cal check  
#7-9

Form 106-I8000