



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Ryan Pomerantz</i>		
Location <i>Minot AFB</i>	Intoxilyzer® Serial Number <i>80-006501</i>	
Gas Lot Number <i>25024080A3</i>	Gas Cylinder Number <i>316</i>	Gas Expiration Date <i>5/12/2028</i>

Check When Done:

- ☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- AEI ☒ 2. Perform an ACA Test (Level 1, Function C).
- 01Oct26 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
- A. *0.080* AC
- B. *0.099* AC
- C. *0.099* AC
- ☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- ☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ☒ 5. Send the Following to the Crime Laboratory:
- A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
- B. ACA Test Record.

Chemical Test Operator Signature <i>Ryan Pomerantz</i>	Date <i>30 Sept 26</i>
Reviewed By (Crime Laboratory Use Only) <i>Anna Ingemansen</i>	Date <i>01 Oct 2026</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006501
Location = MNAF 8164.14.00 09/16
09/30/2026 14:02

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:03
02 Std. Gas	0.080	14:03
03 Room Air	0.000	14:04
04 Std. Gas	0.079	14:04
05 Room Air	0.000	14:05
06 Std. Gas	0.079	14:05
07 Room Air	0.000	14:06

Lot No = 264944
Cyl No = 316
Exp Date = 05/12/2028
County = 64 Oper No. = 134284



Operator Signature
RYAN POMERANZ

Remarks:

Form 106-I8000