

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004193

Location: Statesman Co Corr - upstairs (adult) *JAME*

- Check When Done
- A. Is the warm-up time less than 20 minutes? ✓
 - B. Is a three-pronged grounded outlet used? ✓
 - C. Is the breath tube heated? ✓
 - D. Is the diagnostic check complete? ✓
 - E. Is the time, date, and year correct (re-set if necessary)? ✓
 - F. Print test. (Attach test record.) ✓
 - G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)
Sim SN DR5132 Lot No. 12010 Rep AC 0.020 ✓
 - H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)
Sim SN DR3803 Lot No. 12050 Rep AC 0.300 ✓
 - I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.
(Attach test record.)
Sim SN DR3846 Lot No. Rep AC 0.100 ✓
 - J. RFI check. Run CMS mode. Key radio on first room air.
(Attach test record.) ✓
 - K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.
(Attach test record.)

Lot No. <u>05911080A1</u>	Cyl No. <u>21</u>	Exp. Date <u>02/01/2013</u>
Test 1 <u>0.078</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.078</u>	Test 5 <u>0.079</u>	Test 8 <u>0.079</u>
Test 3 <u>0.078</u>	Test 6 <u>0.078</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 491 psi Reg 500 psi

M. Remarks/Maintenance record: Annual inspection

16 May 2012
Date
17 May 12
Date

Deb Kashur
Field Inspector's Signature
Charles E. Edr
Reviewed by

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:28

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNopqrstuvwxyz!@#%&^&*()_+?

Current Instrument Setup

Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 05911080A1
Standard Cyl #? 21
Standard Expiration? 02/01/2013
Oper No? 777777

Flow Cal. Date: 06/04/2009
Slope 692
Intercept -611138

IR Calibration Date: 05/24/2011
 3um 9um

0th Coef(*100): -19796 -18666
1st Coef(*100): 277921 137858
2nd Coef(*100): -831 322
H2O adj(mg/l*10k): 446 455

***** Printer Test End *****

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Print test

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:33

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:33
02 Std. Sol.	0.021	17:35
03 Room Air	0.000	17:36
04 Std. Sol.	0.021	17:37
05 Room Air	0.000	17:37
06 Std. Sol.	0.021	17:38
07 Room Air	0.000	17:38

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 47

Oper No. = 777777

Deb Kashur

Operator Signature
DEB KASHUR

Remarks:

Low AC 0.020 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:40

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:40
02 Std. Sol.	0.302	17:41
03 Room Air	0.000	17:42
04 Std. Sol.	0.303	17:42
05 Room Air	0.000	17:43
06 Std. Sol.	0.302	17:44
07 Room Air	0.000	17:44

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 47

Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Linearity test ^{DK} ~~at~~ 0.300 AC

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:46

Test	AC	Time
01 Room Air	0.000	17:47
02 *Subject Test	INT*	17:47
03 Room Air	0.000	17:48

*Invalid Test
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 05911080A1
Cyl No = 21
Expiration Date = 02/01/2013
County = 47 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Interferent check 0.100 AC + acetone

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:48

Test	AC	Time
01 Diagnostic	OK	17:48
02 Room Air	RFI*	17:49
03 Room Air	0.000	17:49

*Invalid Test
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT
Sub DOB = 02/01/1992
Sub Sex = Female Weight = 150
Test = DUI Cit = NA
Dr. Lic. = ND/DIS921456
Lot No = 05911080A1
Cyl No = 21
Expiration Date = 02/01/2013
County = 47 Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.



Operator Signature
DEB KASHUR

Remarks:

RFI check

Form 106-I8000

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:29
02 Std. Gas	0.078	17:30
03 Room Air	0.000	17:30
04 Std. Gas	0.078	17:30
05 Room Air	0.000	17:31
06 Std. Gas	0.078	17:31
07 Room Air	0.000	17:32

Lot No = 05911080A1
Cyl No = 21
Exp Date = 02/01/2013
County = 47 Oper No. = 777777

Deb Kasher

Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:51
02 Std. Gas	0.079	17:51
03 Room Air	0.000	17:51
04 Std. Gas	0.079	17:52
05 Room Air	0.000	17:52
06 Std. Gas	0.078	17:53
07 Room Air	0.000	17:53

Lot No = 05911080A1
Cyl No = 21
Exp Date = 02/01/2013
County = 47 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = JAME 8164.13.00 06/09
05/16/2012 17:53

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:54
02 Std. Gas	0.079	17:54
03 Room Air	0.000	17:55
04 Std. Gas	0.079	17:55
05 Room Air	0.000	17:56
06 Std. Gas	0.079	17:56
07 Room Air	0.000	17:56

Lot No = 05911080A1
Cyl No = 21
Exp Date = 02/01/2013
County = 47 Oper No. = 777777

Deb Kashur
Operator Signature
DEB KASHUR

Remarks:

Calibration check 0.080 AC