



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Darcy W Aberle</u>		
Location: <u>Williston LEC</u>		
Intoxilyzer® Serial Number: <u>80-004184</u>		
Gas Lot Number: <u>15410080A1</u>	Gas Cylinder Number: <u>6</u>	Expiration Date: <u>07/01/2012</u>
<input type="checkbox"/> New Cylinder Installed <input checked="" type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%) Test 1: <u>0.080</u> % Test 2: <u>0.079</u> % Test 3: <u>0.080</u> %		
(ATTACH TEST RECORD)		

Darcy W Aberle
 Chemical Test Operator Signature

02/16/2011
 Date

Deb Kashner
 Reviewed By

03 Mar 2011
 Date

80-004184
 P. 102 of 103
 KM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004184
Location = WILL 8164.13.00 06/09
02/16/2011 20:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:36
02 Std. Gas	0.080	20:36
03 Room Air	0.000	20:37
04 Std. Gas	0.079	20:37
05 Room Air	0.000	20:38
06 Std. Gas	0.080	20:38
07 Room Air	0.000	20:39

Lot No = 15410080A1
Cyl No = 6
Exp Date = 07/01/2012
County = 53 Oper No. = 010101


Operator Signature
DARCY W. ABERLE

Remarks:

Form 106-I8000

80-004184
P.103 of 103
VM