



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <u>Robert Kennedy</u>		
Location <u>Devils Lake LEC</u>	Intoxilyzer® Serial Number <u>80-004951</u>	
Gas Lot Number <u>13374080A3</u>	Gas Cylinder Number <u>36</u>	Gas Expiration Date <u>6/5/2016</u>

Check When Done:

☒ 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

☒ 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0. 080 AC

B. 0. 079 AC

C. 0. 080 AC

☒ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

☒ 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

☒ 5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <u>[Signature]</u>	Date <u>7/25/2014</u>
Reviewed By (Crime Laboratory Use Only) <u>Deb Kashur</u>	Date <u>30 July 14</u>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004951
Location = DEVL 8164.13.00 06/09
07/25/2014 11:58

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:59
02 Std. Gas	0.080	11:59
03 Room Air	0.000	12:00
04 Std. Gas	0.079	12:00
05 Room Air	0.000	12:00
06 Std. Gas	0.080	12:01
07 Room Air	0.000	12:01


Lot No = 13314080A3

Cyl No = 36

Exp Date = 06/05/2016

County = 36

Oper No. = 069108



Operator Signature
ROBERT KENNEDY

Remarks:

Form 106-I8000