



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Lawson, Johnny Z</u>		
Location: <u>Bowman SO</u>		
Intoxilyzer® Serial Number: <u>80-004188</u>		
Gas Lot Number: <u>16512080A2</u>	Gas Cylinder Number: <u>31</u>	Expiration Date: <u>08/01/2014</u>
<input type="checkbox"/> New Cylinder Installed <input checked="" type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>077</u> %		
Test 2: 0. <u>078</u> %		
Test 3: 0. <u>077</u> %		
(ATTACH TEST RECORD)		

[Signature]
Chemical Test Operator Signature

08/23/2013
Date

Deb Kashner
Reviewed By

30 Aug. 13
Date


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = BOWM 8164.13.00 06/09
08/23/2013 09:36

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:37
02 Std. Gas	0.077	09:37
03 Room Air	0.000	09:37
04 Std. Gas	0.078	09:38
05 Room Air	0.000	09:38
06 Std. Gas	0.077	09:39
07 Room Air	0.000	09:39

Lot No = 16512080A2
Cyl No = 31
Exp Date = 08/01/2014
County = 06 Oper No. = 020814


Operator Signature
JOHNNY LAWSON

Remarks: