



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Terry Dilloe</i>		
Location <i>CARS</i>	Intoxilyzer® Serial Number <i>90-007094</i>	
Gas Lot Number <i>13518080A6</i>	Gas Cylinder Number <i>019</i>	Gas Expiration Date <i>8-5-2020</i>

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

*ACA test was performed. See pg 2. 17 Jan 19* Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0.079 AC

B. 0.079 AC

C. 0.079 AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <i>Terry Dilloe</i>	Date <i>1-3-19</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Muegge Numms</i>	Date <i>1-17-19</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007094  
Location = CARS      8164.16.00 09/18  
01/03/2019      15:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:41
02 Std. Gas	0.079	15:41
03 Room Air	0.000	15:41
04 Std. Gas	0.079	15:42
05 Room Air	0.000	15:42
06 Std. Gas	0.079	15:43
07 Room Air	0.000	15:43

Lot No = 13518080A6  
Cyl No = 19  
Exp Date = 08/05/2020  
County = 19      Oper No. = 132318

  
Operator Signature  
TERRY DILLON

Remarks:

Form 106-I8000