







**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Wasley Straight</i>		
Location <i>Coop</i>	Intoxilyzer® Serial Number <i>80-004952</i>	
Gas Lot Number <i>24119080A1</i>	Gas Cylinder Number <i>4</i>	Gas Expiration Date <i>11/5/21</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. *0.079* AC  
 B. *0.078* AC  
 C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:  
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
 B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i>	Date <i>7/9/21</i>
Reviewed By (Crime Laboratory Use Only)	Date

NOTE:

EMAILED WES STRAIGHT ON 7.14.20  
 TO CORRECT DATE ON FORM 105-G.

*CEE 7.15.20*