

✓
Dat above
cme

**INTOXILYZER® 5000 G AND 8000
INSTALLATION AND REPAIR CHECKOUT**

Serial No.: 80-004193 Instrument Location: GRAF

Reason for Install/Repair: _____ Install After Lab Repair/Inspection
X Other (specify) NEW INSTALL

Check When Done:

- ✓ A. Surge Protector Installed/Properly Grounded.
- (N/A)✓ B. Telephone Line Connected to Intoxilyzer® KB. NOT AVAILABLE AT TIME OF INSTALL
- ✓ C. Breath Tube Heated.
- ✓ D. Review/Enter Preliminary Data Entry (Level 2, Function E).
- ✓ E. Ethanol Breath Standard: X Connected X Installed.
(Include Form 105-G if New Cylinder Installed.)
- ✓ F. Review/Enter Gas Setup (Level 1, Function S). (Intoxilyzer® 8000 Only.)
- ✓ G. Print Test (Level 1, Function P). (Attach Test Record.)
- ✓ H. ACA Test (Level 1, Function C). (Attach Test Record.)
- ✓ I. RFI Test (Use CMS Mode-First Room Air). (Attach Test Record.)
- J. Repair and/or Maintenance Performed: _____

10/09/09
Date
16 Oct 09
Date

[Signature]
Field Inspector's Signature
Deb Kashur
Reviewed By

80-004193
p. 1 of 4
✓

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = GRAF 8164.12.00 06/09
10/09/2009 10:01

***** Printer Test *****

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=_|
ABCDEFGHIJKLMNPOQRSTUVWXYZ!@#%&^*()_+?

Current Instrument Setup


Data Entry Mode: Enabled
Start Test Sequence: DABACABA
Display Prelim Rslt? Yes
Display Third Digit? Yes
Inhib Printer(Y/N)? No
Display Volume? No
Disable On Memfull? Yes
of Print Copies? 1
Select Std (D/W/I)? Dry
Standard Value? 0.080
Standard Lot #? 659358
Standard Cyl #? 16
Standard Expiration? 09/03/2011
Oper No? 048915

Flow Cal. Date: 06/04/2009
Slope 692
Intercept -611138

IR Calibration Date: 06/03/2009
 3um 9um

0th Coef(*100): -12132 -15777
1st Coef(*100): 261160 129907
2nd Coef(*100): 1958 971
H2O adj(mg/l*10k): 413 335

***** Printer Test End *****



Operator Signature
DOLF H. OLDENBURG

Remarks:

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = GRAF 8164.12.00 06/09
10/09/2009 10:02

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 10:03 |
| 02 Std. Gas | 0.080 | 10:03 |
| 03 Room Air | 0.000 | 10:04 |
| 04 Std. Gas | 0.081 | 10:04 |
| 05 Room Air | 0.000 | 10:04 |
| 06 Std. Gas | 0.080 | 10:05 |
| 07 Room Air | 0.000 | 10:05 |

Lot No = 659358
Cyl No = 16
Exp Date = 09/03/2011
County = 50 Oper No. = 048915

D. Oldenburg 774 NDIA

Operator Signature
DOLF H. OLDENBURG

Remarks:

Form 106-I8000

80-004193
P. 3 of 4 PK

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004193
Location = GRAF 8164.12.00 06/09
10/09/2009 10:07

DRY CAL CHECK

| Test | AC | Time |
|--------------|-------|-------|
| 01 Room Air | 0.000 | 10:08 |
| 02 Reference | RFI* | 10:08 |
| 03 Room Air | 0.000 | 10:08 |

*Invalid Test
Inhibited - RFI

Lot No = 659358

Cyl No = 16

Exp Date = 09/03/2011

County = 50

Oper No. = 048915

D. Oldenburg 774 NDHP.

Operator Signature
DOLF H. OLDENBURG

Remarks:

Form 106-I8000

80-004193
P. 4 of 4 FKL