



**ETHANOL GAS STANDARD CYLINDER REPORT**  
NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Eric Romero</i>		
Location <i>T. Ogil</i>	Intoxilyzer® Serial Number <i>80-007098</i>	
Gas Lot Number <i>02821080A1</i>	Gas Cylinder Number <i>19</i>	Gas Expiration Date <i>3-5-23</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.082 AC
  - B. 0.082 AC
  - C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>[Signature]</i> 1139	Date <i>5/31/22</i>
Reviewed By (Crime Laboratory Use Only) <i>Janelle Pertschiller</i>	Date <i>01 June 22</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-007098  
Location = TIOG      8164.16.00 09/18  
05/31/2022      20:48

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:49
02 Std. Gas	0.082	20:49
03 Room Air	0.000	20:50
04 Std. Gas	0.082	20:50
05 Room Air	0.000	20:51
06 Std. Gas	0.081	20:51
07 Room Air	0.000	20:52

Lot No = 02621080A1  
Cyl No = 19  
Exp Date = 03/05/2023  
County = 53      Oper No. = 132508



Operator Signature  
ERIC ROMERO

Remarks:

Form 106-I8000