

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

| Chemical Test Operator Name (Print) Jeffrey Bauske | | | |
|---|------------------------|---|-----------------------------------|
| Location Devils Lake Law Enforcement Center | | Intoxilyzer® Serial Number 80-007096 | |
| Gas Lot Number 24119080A1 | Gas Cylinder Number 27 | | Gas Expiration Date 11/05/2021 |

Check When Done:

- X 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- X 2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

- A. 0.050 AC
- B. 0.080 AC
- C. 0.650 AC

X 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

- X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- X 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| Chemical Test Operator Signature | Jeffrey E | Bauskes Date 1/28/2020 |
|---|-----------|------------------------|
| Reviewed By (Crime Laboratory Use Only) | alid | Date 1/28/20 |

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

| CMI, Inc. Intoxilyzer | Alcohol Analyzer |
|-------------------------|------------------|
| North Dakota Model 8000 | SN 80-007096 |
| Location = DEVL | 8164.16.00 09/18 |
| 01/28/2020 | 11:59 |

| | | | DRY CAL CHECK | |
|----|------|-----|---------------|-------|
| Te | est | | AC | Time |
| 01 | Room | Air | 0.000 | 12:00 |
| 02 | Std. | Gas | 0.080 | 12:00 |
| 03 | Room | Air | 0.000 | 12:01 |
| 04 | Std. | Gas | 0.080 | 12:01 |
| 05 | Room | Air | 0.000 | 12:01 |
| 06 | Std. | Gas | 0.080 | 12:02 |
| 07 | Room | Air | 0.000 | 12:02 |

-Operator Signature

JEFFREY BAUSKE

Remarks:

Form 106-18000