

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Collins, Travis L.			
Location		Intoxilyzer® Serial Number	
HETT		80-007095	
Gas Lot Number	Gas Cylinder Number	Gas Expiration Date	
24119080A1	25	11/05/2021	

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- Z. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.08 2 AC
 - B. 0.482 AC
 - C. 04 82 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 🗵 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ★ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	Date 1/31/2020
Reviewed By (Crime Laboratory Use Only)	Date 1/31/20

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = HETT
01/31/2020
Alcohol Analyzer
SN 80-007095
8164.16.00 09/18

DRY CAL CHECK

Test			AC	Time	
02 03 04 05 06	Room Std. Room Std. Room Std. Room	Gas Air Gas Air Gas	0.000 0.082 0.000 0.082 0.000 0.082	09:42 09:42 09:42 09:43 09:43 09:44	
				02.44	

Lot No = 24119080A1

Cyl No = 25

Exp Date = 11/05/2021

County = 01

'Oper No. = 020533

Operator Signature TRAVIS COLLINS

Remarks:

Form 106-I8000