



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Collins, Travis		
Location HETT	Intoxilyzer® Serial Number 80-007095	
Gas Lot Number 02621080A1	Gas Cylinder Number 3	Gas Expiration Date 03/05/2023

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. 0.082 AC
 B. 0.082 AC
 C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 B. ACA Test Record.

Chemical Test Operator Signature <i>T Collins</i>	Date 11/04/2021
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ed</i>	Date 11.4.2021


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-007095
Location = HETT 8164.16.00 09/18
11/04/2021 08:23

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:23
02 Std. Gas	0.082	08:24
03 Room Air	0.000	08:24
04 Std. Gas	0.082	08:24
05 Room Air	0.000	08:25
06 Std. Gas	0.081	08:25
07 Room Air	0.000	08:26

Lot No = 02621080A1
Cyl No = 3
Exp Date = 03/05/2023
County = 01 Oper No. = 020533



Operator Signature
TRAVIS COLLINS

Remarks: *ACA - cylinder change*

Form 106-I8000