



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Sarah D Warner</i>		
Location <i>Mott</i>	Intoxilyzer® Serial Number <i>80-007091</i>	
Gas Lot Number <i>26021080A1</i>	Gas Cylinder Number <i>9</i>	Gas Expiration Date <i>10/5/2023</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. *0.079* AC
  - B. *0.080* AC
  - C. *0.079* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Sarah D Warner</i>	Date <i>4/6/2022</i>
Reviewed By (Crime Laboratory Use Only) <i>Janelle Putschler</i>	Date <i>11 May 2022</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-007091  
Location = MOTT                    8164.16.00 09/18  
04/06/2022                            15:35

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	15:36
02 Std. Gas	0.079	15:36
03 Room Air	0.000	15:37
04 Std. Gas	0.080	15:37
05 Room Air	0.000	15:38
06 Std. Gas	0.079	15:38
07 Room Air	0.000	15:39

Lot No = 26021080A1

Cyl No = 9

Exp Date = 10/05/2023

County = 21

Oper No. = 100620



Operator Signature  
SARAH WARNER

Remarks:

Form 106-I8000