

CHAIN OF CUSTODY
BREATH ALCOHOL INSTRUMENT

Instrument: 8000 Serial Number: 80-006673

Agency being assigned instrument: OFFICE OF ATTORNEY GENERAL - CRIME LAB

Date 4-6-2017

Transferred from:

Name (print) CMI, INC.

Signature _____

Agency CMI, INC.

Transferred to:

Name (print) CHARLES E. EDER

Signature Charles E. Eder

Agency O.A.G. - CRIME LAB

Reason for transfer/remarks:

- RETURNED AFTER WARRANTY REPAIR AT CMI
- BATTERY CONNECTOR ADJUSTED/RE-SEATED

Or Shipped Via:

USPS _____ UPS _____ FedEx _____ Other _____

Date Sending Instrument: _____ / _____ / _____ Air Bill No. _____

SERVICE WORK ORDER

1000 E. 9th St., Owensboro KY 42303
 Phone 866-835-0690
 Fax: 270-685-6268

340541

BILL TO:

ATTN:

CITY STATE ZIP

CUSTOMER NO. **585013**

PHONE

SHIP TO: *ND Crime Lab*

ATTN: *Deb Shanaver*

2641 E. Main Ave.

CITY STATE ZIP
Bismarck ND 58501

PHONE
701-328-1619

FAX
701-328-1685

MODEL # **8000** SERIAL # **10-00161673** BILL CODE **3**

EXTRA PARTS RCVD
BT gas box battery, PC

DESCRIPTION OF PROBLEM
Will not run on battery.

Infrared / Fuel Cell
 PRELIMINARY TEST INFORMATION

BATTERY- EXT *N/A* V; INT *N/A* V

external battery connector on the instrument side needed sealed all the way VERIFIED PROBLEM? *Yes*

WORK PERFORMED:
*Re-sealed battery connector.
 Recalibrated & finally tested.*

Miscellaneous
 PRELIMINARY TEST INFORMATION

WORK PERFORMED:
[Signature]

PARTS USED		
QTY	PART NO.	DESCRIPTION

REPAIR	CAL. ADJ.	CALIBRATION	FINAL
<i>mv</i>	<i>ADP</i>	<i>ADP</i>	<i>ADP</i>
½ HOURS	HOURS	HOURS	HOURS
<i>6.0</i>	<i>0.5</i>	<i>1.0</i>	<i>1.0</i>
DATE	DATE	DATE	DATE
<i>Mar 3 2017</i>	<i>Mar 30 17</i>	<i>Mar 30 17</i>	<i>Mar 30 17</i>

SHIPPING STICKER HERE

 SPECIAL SHIPPING INSTRUCTIONS