



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Ray Weber</i>		
Location <i>Hill CTRAIL county S.O.</i>	Intoxilyzer® Serial Number <i>80-006678</i>	
Gas Lot Number <i>03417080A3</i>	Gas Cylinder Number <i>#33</i>	Gas Expiration Date <i>03/05/2019</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. 0. *080* AC
 B. 0. *080* AC
 C. 0. *080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Ray Weber</i>	Date <i>09/12/2017</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Hueque Nunez</i>	Date <i>09/14/2017</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006678
Location = HILL 8164.14.00 09/16
09/12/2017 14:37

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:38
02 Std. Gas	0.080	14:38
03 Room Air	0.000	14:39
04 Std. Gas	0.080	14:39
05 Room Air	0.000	14:40
06 Std. Gas	0.080	14:40
07 Room Air	0.000	14:40

Lot No = 03417080A3
Cyl No = 33
Exp Date = 03/05/2019
County = 49 Oper No. = 058405



Operator Signature
RAY WEBER

Remarks: