



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Kopp, Rebecca</i>		
Location <i>BMDC</i>	<i>Burleigh Norton County Detention Center</i> Instrument Code <i>BMDC</i>	Intoxilyzer® Serial Number <i>50-006674</i>
Gas Lot Number <i>03417080A3</i>	Gas Cylinder Number <i>#35</i>	Gas Expiration Date <i>3-5-19</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 079 AC
 - B. 0. 078 AC
 - C. 0. 078 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Rebecca Kopp</i>	Date <i>9-14-17</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Huegel-Thumme</i>	Date <i>9-20-17</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006674
Location = BMCD 8164.14.00 09/16
09/14/2017 10:03

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:04
02 Std. Gas	0.079	10:04
03 Room Air	0.000	10:05
04 Std. Gas	0.078	10:05
05 Room Air	0.000	10:06
06 Std. Gas	0.078	10:06
07 Room Air	0.000	10:07

Lot No = 03417080A3
Cyl No = 35
Exp Date = 03/05/2019
County = 08 Oper No. = 041005



Operator Signature
REBECCA KOPP

Remarks:

Form 106-I8000