



ETHANOL GAS STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Ervin L. Charette</i>		
Location <i>Rolette County Sheriff</i>		Intoxilyzer® Serial Number <i>80-006511</i>
Gas Lot Number <i>03417086A3</i>	Gas Cylinder Number <i>12</i>	Gas Expiration Date <i>03/05/2019</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.082* AC
 - B. *0.081* AC
 - C. *0.081* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Ervin L. Charette</i>	Date <i>12-17-17</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Inequi-Tummo</i>	Date <i>1-3-18</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006511
Location = ROLL 8164.14.00 09/16
12/17/2017 08:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	08:28
02 Std. Gas	0.082	08:29
03 Room Air	0.000	08:29
04 Std. Gas	0.081	08:29
05 Room Air	0.000	08:30
06 Std. Gas	0.081	08:30
07 Room Air	0.000	08:31


Lot No = 03417080A3

Cyl No = 12

Exp Date = 03/05/2019

County = 40

Oper No. = 130221


Operator Signature
ERVIN CHARETTE

Remarks:

ACA Test

Form 106-I8000