



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) MATTHEW ANDERSON		
Location Richland County Jail	Intoxilyzer® Serial Number 80-006508	
Gas Lot Number 19817080 A4	Gas Cylinder Number 24	Gas Expiration Date 09-05-2019

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.082 AC
 - B. 0.081 AC
 - C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Matthew Anderson</i>	Date 1-31-2018
Reviewed By (Crime Laboratory Use Only) <i>Robatahugue-Norris</i>	Date + 2/5/2018

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006508
Location = WAHP 8164.13.00 06/09
01/31/2018 14:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:51
02 Std. Gas	0.082	14:51
03 Room Air	0.000	14:52
04 Std. Gas	0.081	14:52
05 Room Air	0.000	14:53
06 Std. Gas	0.081	14:53
07 Room Air	0.000	14:53

Lot No = 19817080A4

Cyl No = 24

Exp Date = 09/05/2019

County = 39

Oper No. = 020401



Operator Signature
MATTHEW ANDERSON

Remarks:

Form 106-I8000