



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Allen Kettlyn M</i>		Op# <i>090801</i>
Location <i>Cass Jail</i>		Intoxilyzer® Serial Number <i>80-006668</i>
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>20</i>	Gas Expiration Date <i>2-5-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. 080 AC  
 B. 0. 080 AC  
 C. 0. 079 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>AK 3687</i>	Date <i>7/12/2018</i>
Reviewed By (Crime Laboratory Use Only) <i>Robert [Signature]</i>	Date <i>7/16/2018</i>

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006668  
Location = CASS      8164.14.00 09/16  
07/12/2018      19:27

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	19:27
02 Std. Gas	0.080	19:28
03 Room Air	0.000	19:28
04 Std. Gas	0.080	19:28
05 Room Air	0.000	19:29
06 Std. Gas	0.079	19:29
07 Room Air	0.000	19:30

Lot No = 34917080A3

Cyl No = 20

Exp Date = 02/05/2020

County = 09

Oper No. = 090801

*KA3687*

Operator Signature  
KATHRYN ALLEN

Remarks:

*New gas cylinder instal*

Form 106-I8000