



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Allen Kathryn Marie</i>		OP # <i>090801</i>
Location <i>Cass County Jail</i>		Intoxilyzer® Serial Number <i>80-006668</i>
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>22</i>	Gas Expiration Date <i>08-05-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 080 AC
 - B. 0. 080 AC
 - C. 0. 080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>K Allen #3687</i>	Date <i>8-28-2018</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Muege-Tummo</i>	Date <i>8-28-2018</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006668
Location = CASS 8164.14.00 09/16
08/28/2018 07:15

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	07:15
02 Std. Gas	0.080	07:16
03 Room Air	0.000	07:16
04 Std. Gas	0.080	07:17
05 Room Air	0.000	07:17
06 Std. Gas	0.080	07:17
07 Room Air	0.000	07:18

Lot No = 34917080A3
Cyl No = 22
Exp Date = 02/05/2020
County = 09 Oper No. = 090801


Operator Signature
KATHRYN ALLEN

Remarks: