



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Allen Kathryn M</i>		
Location <i>Class County Jail</i>		Intoxilyzer® Serial Number <i>80-0066668</i>
Gas Lot Number <i>13518080 AL6</i>	Gas Cylinder Number <i>15</i>	Gas Expiration Date <i>8-5-2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.080 AC*
 - B. *0.082 AC*
 - C. *0.082 AC* ^{RAW} _{15 Oct 18} *0.083 AC*
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Kate Zebitz</i>	Date <i>10-14-18</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Huegu. Nunez</i>	Date <i>10-15-18</i>


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006668
Location = CASS 8164.14.00 09/16
10/14/2018 22:54

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	22:54
02 Std. Gas	0.082	22:55
03 Room Air	0.000	22:55
04 Std. Gas	0.082	22:55
05 Room Air	0.000	22:56
06 Std. Gas	0.083	22:56
07 Room Air	0.000	22:57

Lot No = 13518080A6
Cyl No = 15
Exp Date = 08/05/2020
County = 09 Oper No. = 090801



Operator Signature
KATHRYN ALLEN

Remarks:

New canister installed