



ETHANOL GAS STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (10-2013)

Chemical Test Operator Name (Print) <i>Sarah Warner</i>		
Location <i>Mott</i>	Intoxilyzer® Serial Number <i>80-006513</i>	
Gas Lot Number <i>34917080A3</i>	Gas Cylinder Number <i>18</i>	Gas Expiration Date <i>2/5/2020</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0. 080 AC
 - B. 0. 080 AC
 - C. 0. 080 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Sarah Warner</i>	Date <i>11/2/18</i>
Reviewed By (Crime Laboratory Use Only) <i>Roberta Huggins-Rumour</i>	Date <i>11/2/18</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006513
Location = MOTT 8164.13.00 06/09
11/02/2018 13:55

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:56
02 Std. Gas	0.080	13:56
03 Room Air	0.000	13:57
04 Std. Gas	0.080	13:57
05 Room Air	0.000	13:58
06 Std. Gas	0.080	13:58
07 Room Air	0.000	13:59

Lot No = 34917080A3
Cyl No = 18
Exp Date = 02/05/2020
County = 21 Oper No. = 100620



Operator Signature
SARAH WARNER

Remarks: