

## ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print)  FDE2, CHARCES E	STATE TO	ox, ANNUAR	TUSPECTION
WARD CONNTY S.O.	(JAIL)	Intoxilyzer® Serial Nu	1mber 80 - 006677
Gas Lot Number 344/8080A2	Gas Cylinder Number	48	Gas Expiration Date 2/5/2021

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C). Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0.**282** AC B. 0.**282** AC C. 0.**282** AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

Send the Following to the Crime Laboratory:

A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature	Date 5.6.19
Reviewed By (Crime Laboratory Use Only)	Date 5.6.19

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = WARD
05/06/2019

Alcohol Analyzer SN 80-006677 8164.14.00 09/16 14:00

## DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:01
02 Std. Gas	0.082	14:01
03 Room Air	0.000	14:02
04 Std. Gas	0.082	14:02
05 Room Air	0.000	14:03
06 Std. Gas	0.082	14:03
07 Room Air	0.000	14:04

Lot No = 34418080A2

Cyl No = 48

Exp Date = 02/05/2021

County = 51

Oper No. = 666666

Operator Signature CHARLES EDER

Remarks:

Form 106-I8000