



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) Darcy Aberle		
Location Williams County Corrections		Intoxilyzer® Serial Number 80-006509
Gas Lot Number 24119080A1	Gas Cylinder Number 37	Gas Expiration Date 11/05/21

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. **0.080**AC
 B. **0.080**AC
 C. **0.080**AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 B. ACA Test Record.

Chemical Test Operator Signature Darcy Aberle	Date 11/14/19
Reviewed By (Crime Laboratory Use Only) Charles E. Ehn	Date 11.15.2019


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006509
Location = WILL 8164.14.00 09/16
11/14/2019 16:34

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	16:35
02 Std. Gas	0.080	16:35
03 Room Air	0.000	16:36
04 Std. Gas	0.080	16:36
05 Room Air	0.000	16:36
06 Std. Gas	0.080	16:37
07 Room Air	0.000	16:37

Lot No = 24119080A1
Cyl No = 37
Exp Date = 11/05/2021
County = 53 Oper No. = 010101


Operator Signature
DARCY ABERLE

Remarks:

Form 106-I8000