



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <b>MATTHEW ANDERSON</b>		
Location <b>RICHLAND COUNTY JAIL</b>		Intoxilyzer® Serial Number <b>80 - 006508</b>
Gas Lot Number <b>34418080 A 2</b>	Gas Cylinder Number <b>8</b>	Gas Expiration Date <b>02-05-2021</b>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. **080** AC  
 B. 0. **081** AC  
 C. 0. **081** AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:  
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).  
 B. ACA Test Record.

Chemical Test Operator Signature <i>Matthew Anderson</i>	Date <b>5-28-20</b>
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ehl</i>	Date <b>6/4/2020</b>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-006508  
Location = WAHP                8164.13.00 06/09  
05/28/2020                        14:40

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:40
02 Std. Gas	0.080	14:41
03 Room Air	0.000	14:41
04 Std. Gas	0.081	14:42
05 Room Air	0.000	14:42
06 Std. Gas	0.081	14:43
07 Room Air	0.000	14:43

Lot No = 34418080A2  
Cyl No = 8  
Exp Date = 02/05/2021  
County = 39                        Oper No. = 020401



Operator Signature  
MATTHEW ANDERSON

Remarks:

Form 106-I8000