



**ETHANOL GAS STANDARD CYLINDER REPORT**  
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM  
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) O'Brien, Matthew, Michael		
Location Oakes Police Department		Intoxilyzer® Serial Number 80-006680
Gas Lot Number 05620080A1	Gas Cylinder Number 32	Gas Expiration Date 04/05/2022

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
  - A. 0.083 AC
  - B. 0.082 AC
  - C. 0.083 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Matthew O'Brien</i>	Date 02/02/2021
Reviewed By (Crime Laboratory Use Only) <i>Charles E. Ed</i>	Date 2/4/21


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-006680  
Location = OKPD      8164.14.00 09/16  
02/02/2021      10:12

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	10:13
02 Std. Gas	0.083	10:13
03 Room Air	0.000	10:14
04 Std. Gas	0.082	10:14
05 Room Air	0.000	10:15
06 Std. Gas	0.083	10:15
07 Room Air	0.000	10:16

Lot No = 05620080A1  
Cyl No = 32  
Exp Date = 04/05/2022  
County = 11      Oper No. = 130690

  
\_\_\_\_\_  
Operator Signature  
MATTHEW O'BRIEN

Remarks:

Form 106-I8000