



ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) MARK PEARSON		
Location EMMONS COUNTY	Intoxilyzer® Serial Number 80-000089	
Gas Lot Number 05620080A1	Gas Cylinder Number 29	Gas Expiration Date 4/5/2022

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. **0.079** AC
 B. **0.079** AC
 C. **0.079** AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature	Date 2/2/2021
Reviewed By (Crime Laboratory Use Only)	Date 2.8.21

Intoxilyzer Test Record and Checklist
 CORE Crime Lab. Div., Bismarck, ND 58501
 1, Inc. Intoxilyzer Alcohol Analyz.
 North Dakota Model 8000 SN 80-006689
 Calibration = LINT 8164.14.00 09/2021
 02/2021

DRY CAL CHECK

AC	Time
Room Air	0.000
Std. Gas	0.079
Room Air	0.000
Std. Gas	0.079
Room Air	0.000
Std. Gas	0.079
Room Air	0.000

No = 05620080A1
 No = 29
 Date = 04/05/2022
 Qty = 15

Operator Signature
 MARK PEARSON

marks:

Form 106-18000




ETHANOL GAS STANDARD CYLINDER REPORT
 NORTH DAKOTA OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM
 SFN 59282 (08/2018)

Chemical Test Operator Name (Print) MARK PEARSON		
Location EMMONS COUNTY	Intoxilyzer® Serial Number 80-000689	
Gas Lot Number 05620080A1	Gas Cylinder Number 29	Gas Expiration Date 4/5/2022

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 A. **0.079** AC
 B. **0.079** AC
 C. **0.079** AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 B. ACA Test Record.

Chemical Test Operator Signature 	Date 2/2/2021
Reviewed By (Crime Laboratory Use Only)	Date

Intoxilyzer Test Record and Checklist
 ND006 Crime Lab. Div., Bismarck, ND 58501
 CMI, Inc. Intoxilyzer Alcohol Analyzer
 North Dakota Model 8000 SN 80-006689
 Location = LINT 8164.14.00 09/16
 02/02/2021 11:29

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	11:30
02 Std. Gas	0.079	11:30
03 Room Air	0.000	11:31
04 Std. Gas	0.079	11:31
05 Room Air	0.000	11:32
06 Std. Gas	0.079	11:32
07 Room Air	0.000	11:33

Lot No = 05620080A1
 Cyl No = 29
 Exp Date = 04/05/2022
 County = 15 Oper No. = 130117

Operator Signature
 MARK PEARSON

Remarks:

Form 106-18000

FORM 105-G

CEE 2/4/21 OPERATOR NEEDS TO SIGN