

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

Chemical Test Operator Name (Print)	UNER SMITH	132660
Location GRAFTON LEC	Intoxilyzer® Seri	ial Number 80 - 006686
Gas Lot Number 02621080A Ga	as Cylinder Number	Gas Expiration Date 03/05/2023

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.023 AC
 - B. 0083 AC
 - C. 0.08 3 AC
- Z 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

	duelso	Date 02/19/2021
Reviewed By (Crime Laboratory Use Only)	Charles Esta	Date 02.19.2021

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer North Dakota Model 8000 Location = GRAF Alcohol Analyzer SN 80-006686 8164.14.00 09/16 12:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	12:22
02 Std. Gas	0.083	12:22
03 Room Air	0.000	12:23
04 Std. Gas	0.083	12:23
05 Room Air	0.000	12:23
06 Std. Gas	0.083	12:24
07 Room Air	0.000	12:24

Lot No = 02621080A1

Cyl No = 6

02/19/2021

Exp Date = 03/05/2023

County = 50

Oper No. = 132660

Operator Signature
TANNER SMITH

Remarks:

ACA Test

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006686
Location = GRAF 8164.14.00 09/16
02/19/2021 12:16

******* Printer Test ********

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$\%*()_+?

abcdefghijklmnopqrstuvwxyz1234567890-=|
ABCDEFGHIJKLMNOPQRSTUVWXYZ!@#\$%^&*()_+?

Current Instrument Setup Enabled Data Entry Mode: Start Test Sequence: DABACABA Display Prelim Rslt? Yes Yes Display Third Digit? Inhib Printer(Y/N)? No No Display Volume? Disable On Memfull? Yes 1 # of Print Copies? Select Std (D/W/I)? Dry 0.080 Standard Value? 02621080A1 Standard Lot #? Standard Cyl #? 03/05/2023 Standard Expiration? 132660 Oper No?

Flow Cal. Date: 09/07/2016 Slope 770 Intercept -749410

IR Calibration Date: 09/06/2016

3um 9um

Oth Coef(*100): -14519 -24061

1st Coef(*100): 252211 131592

2nd Coef(*100): 3527 1490

H2O adj(mg/l*10k): 522 594

****** Printer Test End *******

Operator Signature
TANNER SMITH

Remarks:

Form 106-I8000

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = GRAF
02/19/2021

Alcohol Analyzer SN 80-006686 8164.14.00 09/16 12:28

DRY CAL CHECK

Test	AC	Time
01 Room Air	RFI*	12:28
02 Room Air	0.000	12:28

*Invalid Test Inhibited - RFI

Lot No = 02621080A1

Cyl No = 6

Exp Date = 03/05/2023

County = 50

Oper No. = 132660

Operator Signature TANNER SMITH

Remarks:

RFI Test

Form 106-I8000