

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

CHRISTOPHER L. ALMLIE			
Location Intoxilyzer® Serial Number SO-CO691			
Gas Lot Number Gas Cylinder Number Gas Expiration O3/05/	Pate /2023		
Check When Done:			
1 Scan/Enter Gas Cylinder Information (Level 1, Function S).			
☑ 2. Perform an ACA Test (Level 1, Function C).			
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).			
A. 0.₫¼೭ AC			
B. 0.0% AC			
C. 0.Q82 AC	to ili mano fon Haa		
☑ 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.			
4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.			
√ 5. Send the Following to the Crime Laboratory:			
A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).B. ACA Test Record.			
Chemical Test Operator Signature	Date 4/10/2021		
Beviewed By (Crime Laboratory Use Only)	Date 4.23 2021		

FORM 105-G

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = KMPD
04/10/2021

SN 80-006691 8164.14.00 09/16 20:21

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:22
02 Std. Gas	0.082	20:22
03 Room Air	0.000	20:23
04 Std. Gas	0.081	20:23
05 Room Air	0.000	20:24
06 Std. Gas	0.082	20:24
07 Room Air	0.000	20:24

Lot No = 02621080A1

Cyl No = 030

Exp Date = 03/05/2023

County = 51

Oper No. = 132664

Operator Signature CHRISTOPHER ALMLIE

Remarks:

Form 106-I8000