

ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM SFN 59282 (08/2018)

| Chemical Test Operator Name (Print) Jeffrey A. Bauske | | | |
|---|---------------------------|---|-----------------------------------|
| Location DEVL | | Intoxilyzer® Serial Number 80-006671 | |
| Gas Lot Number 26021080A1 | Gas Cylinder Number 19 | | Gas Expiration Date 10/05/2023 |

Check When Done:

- 2. Perform an ACA Test (Level 1, Function C).
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. 0.08 1 AC
 - B. 0.681 AC
 - C. 0.081 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- X 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- ∑ 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

| Chemical Test Operator Signature Ostary Bauska | Date 01/22/2022 |
|---|-----------------|
| Beviewed By (Crime Laboratory Use Only) | Date 1/24/2022 |

Intoxilyzer Test Record and Checklist NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer
North Dakota Model 8000
Location = DEVL

01/22/2022

Alcohol Analyzer SN 80-006671 8164.14.00 09/16 16:10

DRY CAL CHECK

| Test | AC | Time |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 16:11 |
| 02 Std. Gas | 0.081 | 16:11 |
| 03 Room Air | 0.000 | 16:12 |
| 04 Std. Gas | 0.081 | 16:12 |
| 05 Room Air | 0.000 | 16:13 |
| 06 Std. Gas | 0.081 | 16:13 |
| 07 Room Air | 0.000 | 16:14 |

Lot No = 26021080A1

Cyl No = 19

Exp Date = 10/05/2023

County = 36

Oper No. = 120501

Operator Signature JEFFREY BAUSKE

Remarks:

DEA TEST / NEW GAS CYCINDER INSTALL

Form 106-I8000