



ETHANOL GAS STANDARD CYLINDER REPORT

NORTH DAKOTA OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION-TOXICOLOGY SECTION/BREATH ALCOHOL PROGRAM

SFN 59282 (08/2018)

Chemical Test Operator Name (Print) <i>Shayla Rodgers, 4477</i>		
Location <i>Mountrail Co. Booking</i>	Intoxilyzer® Serial Number <i>80-006676 26021080A1 SK</i>	
Gas Lot Number <i>211543</i>	Gas Cylinder Number <i># 016</i>	Gas Expiration Date <i>10/5/2023</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).
Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).
 - A. *0.081* AC
 - B. *0.030* AC
 - C. *0.080* AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
 - A. Ethanol Gas Standard Cylinder Report (SFN59282, Form 105-G).
 - B. ACA Test Record.

Chemical Test Operator Signature <i>Shayla Rodgers 4477</i>	Date <i>6/28/22</i>
Reviewed By (Crime Laboratory Use Only) <i>Janelle Pertschiller</i>	Date <i>7/11/2022</i>

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-006676
Location = MOSO 8164.14.00 09/16
06/28/2022 18:09

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	18:10
02 Std. Gas	0.081	18:10
03 Room Air	0.000	18:10
04 Std. Gas	0.080	18:11
05 Room Air	0.000	18:11
06 Std. Gas	0.080	18:12
07 Room Air	0.000	18:12


Lot No = 26021080A1

Cyl No = 016

Exp Date = 10/05/2023

County = 31

Oper No. = 132957


Operator Signature
SHAYLA RODGERS

Remarks:

Form 106-I8000