



# ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Frederick E. Marquardt</u>		
Location: <u>Grattan LEC</u>		
Intoxilyzer® Serial Number: <u>80-003072</u>		
Gas Lot Number: <u>16</u>	Gas Cylinder Number: <u>659358</u>	Expiration Date: <u>09-03-11</u>
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>83</u> %		
Test 2: 0. <u>82</u> %		
Test 3: 0. <u>82</u> %		
(ATTACH TEST RECORD)		

Frederick E. Marquardt  
Chemical Test Operator Signature

01-13-11  
Date

Deb Kashner  
Reviewed By

24 Jan 2010 <sup>2011 BK</sup>  
Date

80-003072  
P. 100 of 101  
YR


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003072  
Location = GRAF      8164.13.00 06/09  
01/13/2011      02:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	02:51
02 Std. Gas	0.083	02:51
03 Room Air	0.000	02:52
04 Std. Gas	0.082	02:52
05 Room Air	0.000	02:53
06 Std. Gas	0.082	02:53
07 Room Air	0.000	02:53

Lot No = 659358  
Cyl No = 16  
Exp Date = 09/03/2011  
County = 50      Oper No. = 059006

  
Operator Signature  
FREDRICK E. MARQUARDT

Remarks:

*ACA test*

Form 106-I8000

80-003072  
p. 101 of 101  
FM