



ETHANOL BREATH STANDARD CYLINDER REPORT
 OFFICE OF ATTORNEY GENERAL
 CRIME LABORATORY DIVISION
 Toxicology Section/Breath Alcohol Program
 SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Richard D. Johnson</u>		
Location: <u>WASH McLean LEC</u>		
Intoxilyzer® Serial Number: <u>80-004196</u>		
Gas Lot Number: <u>659357</u>	Gas Cylinder Number: <u>47</u>	Expiration Date: <u>9-3-11</u>
<input type="checkbox"/> New Cylinder Installed <input checked="" type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%) Test 1: <u>0.081</u> % Test 2: <u>0.081</u> % Test 3: <u>0.081</u> %		
(ATTACH TEST RECORD)		

Richard D. Johnson
 Chemical Test Operator Signature

1-25-2011
 Date

Deb Kashur
 Reviewed By

28-Jan 2011
 Date

Can you please check the time on intox.
 it is off.

Thanks.

Richard Johnson

80-004196
 P. 99 of 100
 ym

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004196
Location = WASH 8164.13.00 06/09
01/25/2011 13:50

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:51
02 Std. Gas	0.081	13:51
03 Room Air	0.000	13:52
04 Std. Gas	0.081	13:52
05 Room Air	0.000	13:53
06 Std. Gas	0.081	13:53
07 Room Air	0.000	13:53


Lot No = 659357

Cyl No = 47

Exp Date = 09/03/2011

County = 28

Oper No. = 118707



Operator Signature
RICHARD D. JOHNSON

Remarks:

Form 106-I8000

80-004196
P. 100 of 100
VAK