



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>John Stotler</u>		
Location: <u>GFAF</u>		
Intoxilyzer® Serial Number: <u>80-004201</u>		
Gas Lot Number: <u>16512080A2</u>	Gas Cylinder Number: <u>29</u>	Expiration Date: <u>8/1/2014</u>
<input checked="" type="checkbox"/> New Cylinder Installed	<input type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>081</u> %		
Test 2: 0. <u>081</u> %		
Test 3: 0. <u>081</u> %		
(ATTACH TEST RECORD)		

John Stotler
Chemical Test Operator Signature

1/31/2013
Date

Deb Keshner
Reviewed By

25 Feb 2013
Date


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004201
Location = GFAP 8164.13.00 06/09
01/31/2013 20:56

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	20:57
02 Std. Gas	0.081	20:57
03 Room Air	0.000	20:58
04 Std. Gas	0.081	20:58
05 Room Air	0.000	20:59
06 Std. Gas	0.081	20:59
07 Room Air	0.000	21:00

Lot No = 16512080A2
Cyl No = 29
Exp Date = 08/01/2014
County = 65 Oper No. = 031121



Operator Signature
JOHN STOTLER

Remarks: