



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Richard D. Johnson</u>		
Location: <u>McLean LEC WASH</u>		
Intoxilyzer® Serial Number: <u>80-004196</u>		
Gas Lot Number: <u>659357</u>	Gas Cylinder Number: <u>47</u>	Expiration Date: <u>9/30/11</u> ^{DK}
<input type="checkbox"/> New Cylinder Installed <input checked="" type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>082</u> %		
Test 2: 0. <u>082</u> %		
Test 3: 0. <u>081</u> %		
(ATTACH TEST RECORD)		

Richard D. Johnson
Chemical Test Operator Signature

3-1-11
Date

Deb Kashur
Reviewed By

03 Mar 2011
Date

80-004196
P. 101 of 102
FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004196
Location = WASH 8164.13.00 06/09
03/01/2011 13:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:30
02 Std. Gas	0.082	13:31
03 Room Air	0.000	13:31
04 Std. Gas	0.082	13:32
05 Room Air	0.000	13:32
06 Std. Gas	0.081	13:33
07 Room Air	0.000	13:33

Lot No = 659357
Cyl No = 47
Exp Date = 09/03/2011
County = 28 Oper No. = 118707



Operator Signature
RICHARD D. JOHNSON

Remarks:

Form 106-I8000

80-004196
P. 102 of 102
JFM