



**INTOXILYZER® 8000 INSTALLATION AND REPAIR CHECKOUT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59281 (03-2014)

|   |  |
|---|--|
| Serial Number<br><b>80-005938</b>   | Instrument Location<br><b>MOUNTAIN COUNTY JAIL</b> |
| Reason for Install/Repair<br><input checked="" type="checkbox"/> Install After Receiving From Crime Laboratory <input type="checkbox"/> Install After Location Change<br><input type="checkbox"/> Other (Specify) _____ |  |

Check When Done:

- 1. Surge Protector Installed/Properly Grounded.
- 2. Telephone Line Connected to Intoxilyzer® 8000. Available
- 3. Breath Tube Heated.
- 4. Enter Preliminary Data (i.e. Date, Time, DST (Y), and Location; Level 2, Function E).
- 5. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 6. Run Tests:
  - A. Print Test (Level 1, Function P).
  - B. ACA Test (Level 1, Function C).
  - C. RFI Test (CMS Mode or Level 1, Function B or C; Key Radio During Test).
- 7. Repair and/or Maintenance Performed (if Any): \_\_\_\_\_
- 8. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 9. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 10. Send the Following to the Crime Laboratory: Completed Intoxilyzer®8000 Installation and Repair Checkout (SFN59281, Form 104-G), Print Test, ACA Test, and RFI Test.

|  |                            |
|--|----------------------------|
| Field Inspector Signature<br><i>Tom B. Duroe 1235</i>          | Date<br><i>2-3-16 ✓</i>    |
| Reviewed By (Crime Laboratory Use Only)<br><i>Deb Shanaver</i> | Date<br><i>08 Mar 2016</i> |

This record on file at the Office of Attorney General, Crime Laboratory Division, in the County of Burleigh, North Dakota, is certified to be a true and correct copy of the documents received.

Certified By: *Charles E. Edin* Certified Date: *17 Nov. 2016*



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005938  
Location = STAN      8164.13.00 06/09  
03/03/2016      19:59

DRY CAL CHECK

| Test        | AC    | Time  |
|-------------|-------|-------|
| 01 Room Air | 0.000 | 20:00 |
| 02 Std. Gas | 0.081 | 20:00 |
| 03 Room Air | 0.000 | 20:01 |
| 04 Std. Gas | 0.080 | 20:01 |
| 05 Room Air | 0.000 | 20:02 |
| 06 Std. Gas | 0.080 | 20:02 |
| 07 Room Air | 0.000 | 20:03 |

Lot No = 13314080A2  
Cyl No = 054  
Exp Date = 06/05/2016  
County = 31      Oper No. = 031123

  
Operator Signature  
JOHN TRONSET

Remarks:

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-005938  
Location = STAN      8164.13.00 06/09  
03/03/2016      20:04

DRY CAL CHECK

| Test         | AC    | Time  |
|--------------|-------|-------|
| 01 Room Air  | 0.000 | 20:05 |
| 02 Reference | RFI*  | 20:05 |
| 03 Room Air  | 0.000 | 20:05 |

\*Invalid Test  
Inhibited - RFI

Lot No = 13314080A2  
Cyl No = 054  
Exp Date = 06/05/2016  
County = 31      Oper No. = 031123

*John Tronset*  
Operator Signature  
JOHN TRONSET

TEST

Remarks: