



# ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL

CRIME LABORATORY DIVISION

Toxicology Section/Breath Alcohol Program

SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Joseph Knowski</i>		
Location <i>Devils lake LEC</i>	Intoxilyzer® Serial Number <i>80-003061</i>	
Gas Lot Number <i>16512080A2</i>	Gas Cylinder Number <i>08</i>	Gas Expiration Date <i>8-1-2014</i>

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).

A. 0. 081 AC

B. 0. 081 AC

C. 0. 081 AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <i>Joseph Knowski</i>	Date <i>3-8-2013</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>18 Mar 13</i>


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003061  
Location = DEVL      8164.13.00 06/09  
03/08/2013      14:19

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:20
02 Std. Gas	0.081	14:20
03 Room Air	0.000	14:21
04 Std. Gas	0.081	14:21
05 Room Air	0.000	14:22
06 Std. Gas	0.081	14:22
07 Room Air	0.000	14:23

Lot No = 16512080A2  
Cyl No = 8  
Exp Date = 08/01/2014  
County = 36      Oper No. = 109208

  
Operator Signature  
JOSEPH KNOWSKI

Remarks: