

03/11/2013 11:43 ND CRIME LAB

(FAX) 701 328 6185

P.002/003



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Submittal Form (Rev. 10/1/07)
SN: 15228 (3/2011)

Operator Name (Print)	Johnny R. Lawson		
Location	Bowman Co SD	Analyzer & Serial Number	80-004788
Site Log # (e.g. 231M 080A1)	Gas Cylinder Number	Gas Regulation Date	10-01-2013
	69		

Check When Done

- 1. Scan/Enter Gas Cylinder Information (Serial #, Function B)
- 2. Perform an ACA Test (Event I, Function C)
Write in the ACA Test Results (Reported to 3 Digits Ex. 01571)
 - A. 0.679AC
 - B. 0.056AC
 - C. 0.030AC
- 3. Complete the Top Portion of the Worksheet @ Record (SPR2005A, Form 170-6) and place it by the Analyzer @ for file.
- 4. File Previous Worksheet @ Record (SPR2005A, Form 170-6) in the Analyzer @ Location at the Agency
- 5. Send the Worksheet to the Chief Laboratory
 - A. Ethanol Breath Standard Cylinder Report (SPR2005A, Form 100-6)
 - B. ACA Test Report

Chemical Test Operator Signature	Date
<i>[Signature]</i>	03-08-2013
Reviewed by Crime Laboratory Supervisor	Date
<i>[Signature]</i>	18 Mar 13

Form 18-6



ETHANOL BREATH STANDARD CYLINDER REPORT
OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (3-2011)

Chemical Test Operator Name (Print) <u>Johnny Z. Lawson</u>		
Location <u>Bowman Co SD</u>	Intoxilyzer® Serial Number <u>80-004188</u>	
Gas Lot Number <u>83911080A1</u>	Gas Cylinder Number <u>67</u>	Gas Expiration Date <u>10-01-2013</u>

Check When Done:

1. Scan/Enter Gas Cylinder Information (Level 1, Function S).

2. Perform an ACA Test (Level 1, Function C).

Write in the ACA Test Results (Reported to 3 Digits, Ex: 0.081).

A. 0. ___ AC

B. 0. ___ AC

C. 0. ___ AC

3. Complete the Top Portion of the Intoxilyzer® Record (SFNS0496, Form 120-G) and Place it by the Intoxilyzer® for Use.

4. File Previous Intoxilyzer® Record (SFNS0496, Form 120-G) at the Intoxilyzer® Location at the Agency.

5. Send the Following to the Crime Laboratory:

A. Ethanol Breath Standard Cylinder Report (SFNS9282, Form 105-G).

B. ACA Test Record.

Chemical Test Operator Signature <u>[Signature]</u>	Date <u>03-08-2013</u>
Reviewed By (Crime Laboratory Use Only)	Date

Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004188
Location = BOWM 8164.13.00 06/09
03/08/2013 09:20

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:21
02 Std. Gas	0.079	09:21
03 Room Air	0.000	09:22
04 Std. Gas	0.080	09:22
05 Room Air	0.000	09:23
06 Std. Gas	0.080	09:23
07 Room Air	0.000	09:23

Lot No = 23411080A1
Cyl No = 69
Exp Date = 10/01/2013
County = 06

Oper No. = 020814



Operator Signature
JOHNNY LAWSON

Remarks: