



**ETHANOL BREATH STANDARD CYLINDER REPORT**  
OFFICE OF ATTORNEY GENERAL  
CRIME LABORATORY DIVISION  
Toxicology Section/Breath Alcohol Program  
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name: <u>Will Pitt</u>		
Location: <u>GFAF</u>		
Intoxilyzer® Serial Number: <u>80-002667</u>		
Gas Lot Number: <u>05911080A</u>	Gas Cylinder Number: <u>10</u>	Expiration Date: <u>02/01/2013</u>
<input checked="" type="checkbox"/> New Cylinder Installed	<input type="checkbox"/> Cylinder Not Changed	
Instrument Test Results (Report to 3 Digits; Example, 0.081%)		
Test 1: 0. <u>078</u> %		
Test 2: 0. <u>078</u> %		
Test 3: 0. <u>078</u> %		
(ATTACH TEST RECORD)		

Will Pitt  
Chemical Test Operator Signature

21 Mar 11  
Date

Deb Kashur  
Reviewed By

14 Apr 2011  
Date

80-002667  
p. 90 of 91  
RM


Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-002667  
Location = GFAP      8164.13.00 06/09  
03/21/2011      14:38

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	14:39
02 Std. Gas	0.078	14:39
03 Room Air	0.000	14:40
04 Std. Gas	0.078	14:40
05 Room Air	0.000	14:41
06 Std. Gas	0.078	14:41
07 Room Air	0.000	14:42

Lot No = 05911080A1  
Cyl No = 10  
Exp Date = 02/01/2013  
County = 18      Oper No. = 120011



Operator Signature  
WILLIAM R. PETTY

Remarks: Re-install 2nd ACA

Form 106-I8000

80-002667  
P. 9/1/91  
FM