



ETHANOL BREATH STANDARD CYLINDER REPORT

OFFICE OF ATTORNEY GENERAL
CRIME LABORATORY DIVISION
Toxicology Section/Breath Alcohol Program
SFN 59282 (1-2010)

Only Chemical Test Operators may change the Ethanol Breath Standard and complete this form.

Chemical Test Operator Name:	David Feland		
Location:	Kidder S.O.		
Intoxilyzer® Serial Number:	80-004186		
Gas Lot Number:	659357	Gas Cylinder Number:	137 37 ^{OK 14 Apr 11}
		Expiration Date:	09/03/11
<input type="checkbox"/> New Cylinder Installed	<input checked="" type="checkbox"/> Cylinder Not Changed		
Instrument Test Results (Report to 3 Digits; Example, 0.081%)			
Test 1:	0.081	%	
Test 2:	0.080	%	
Test 3:	0.081	%	
			(ATTACH TEST RECORD)

David Feland
Chemical Test Operator Signature

30 Mar 11
Date

Deb Kashner
Reviewed By

14 Apr 11
Date

80-004186
P. 101 of 103
FM


Intoxilyzer Test Record and Checklist
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer Alcohol Analyzer
North Dakota Model 8000 SN 80-004186
Location = STEE 8164.13.00 06/09
03/30/2011 13:25

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	13:26
02 Std. Gas	0.081	13:26
03 Room Air	0.000	13:27
04 Std. Gas	0.080	13:27
05 Room Air	0.000	13:28
06 Std. Gas	0.081	13:28
07 Room Air	0.000	13:29

Lot No = 659357
Cyl No = 37
Exp Date = 09/03/2011
County = 22 Oper No. = 119404



Operator Signature
DAVID FELAND

Remarks:

Form 106-I8000

80-004186
p. 102 of 103
FM

25 April 2011

I called David Feland regarding an error he has been making on 8000 ACA paperwork. He has been writing the gas cylinder number as 137 when in fact it is 37. The correct information is in the Intox, but he writes it incorrectly on the Form 105-G.

Debs Kashner

80-004186
p. 103 of 103
AM