

INTOXILYZER 8000 INSPECTION

Serial No.: 80-004189

Location: TOXL

- |                                                                                                                                                                             | Check When Done |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------|
| A. Is the warm-up time less than 20 minutes?                                                                                                                                | <u>✓</u>        |
| B. Is a three-pronged grounded outlet used?                                                                                                                                 | <u>✓</u>        |
| C. Is the breath tube heated?                                                                                                                                               | <u>✓</u>        |
| D. Is the diagnostic check complete?                                                                                                                                        | <u>✓</u>        |
| E. Is the time, date, and year correct (re-set if necessary)?                                                                                                               | <u>✓</u>        |
| F. Print test. (Attach test record.)                                                                                                                                        | <u>✓</u>        |
| G. Low AC. Use <0.03 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5132</u> Lot No. <u>    </u> Rep AC <u>0.021</u>                                                  | <u>✓</u>        |
| H. Linearity Test. Use ≥ 0.25 AC in ACA mode. (Attach test record.)<br>Sim SN <u>DR5134</u> Lot No. <u>    </u> Rep AC <u>0.294</u>                                         | <u>✓</u>        |
| I. Interferent check. Use 0.05% acetone plus 0.10 AC ethanol in ABA mode.<br>(Attach test record.)<br>Sim SN <u>DR3846</u> Lot No. <u>    </u> Rep AC <u>0.10 + acetone</u> | <u>✓</u>        |
| J. RFI check. Run CMS mode. Key radio on first room air.<br>(Attach test record.)                                                                                           | <u>✓</u>        |
| K. Calibration Check. Use Ethanol Breath Standard cylinder. Do three sets in ACA mode.<br>(Attach test record.)                                                             |                 |

Lot No. <u>05911080A1</u>	Cyl No. <u>17</u>	Exp. Date <u>02/01/2013</u>
Test 1 <u>0.078</u>	Test 4 <u>0.079</u>	Test 7 <u>0.079</u>
Test 2 <u>0.079</u>	Test 5 <u>0.080</u>	Test 8 <u>0.080</u>
Test 3 <u>0.080</u>	Test 6 <u>0.079</u>	Test 9 <u>0.079</u>
Average <u>0.079</u>		

L. Tank pressure, Level 3 Func D Sub G P: 639 psi Reg 650 psi

M. Remarks/Maintenance record: Annual inspection

01 Apr 2011  
Date

Deb Kashur  
Field Inspector's Signature

1 Apr 2011  
Date

Mary Beam  
Reviewed by

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000            SN 80-004189  
Location = TOXL                    8164.13.00 06/09  
04/01/2011                            09:26

\*\*\*\*\* Printer Test \*\*\*\*\*

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

abcdefghijklmnopqrstuvwxyz1234567890-=\_|  
ABCDEFGHIJKLMNopqrstuvwxyz!@#\$\$%^&\*()\_+?

Current Instrument Setup

Data Entry Mode:            Enabled  
Start Test Sequence:        DABACABA  
Display Prelim Rslt?        Yes  
Display Third Digit?        Yes  
Inhib Printer(Y/N)?        No  
Display Volume?            No  
Disable On Memfull?        Yes  
# of Print Copies?          1  
Select Std (D/W/I)?        Dry  
Standard Value?            0.080  
Standard Lot #?            05911080A1  
Standard Cyl #?            17  
Standard Expiration?       02/01/2013  
Oper No?                    777777

Flow Cal. Date:            02/02/2011  
Slope                        628  
Intercept                   -452608

IR Calibration Date:       02/02/2011  
                                    3um                            9um

-----  
0th Coef(\*100):            -22075                        -26801  
1st Coef(\*100):            257433                        132480  
2nd Coef(\*100):            2132                            1061  
H2O adj(mg/l\*10k):        437                            430

\*\*\*\*\* Printer Test End \*\*\*\*\*

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

Print test

Form 106-I8000

80-004189  
P. 113 of 120  
FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = TOXL      8164.13.00 06/09  
04/01/2011      09:54

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:55
02 Std. Sol.	0.019	09:55
03 Room Air	0.000	09:56
04 Std. Sol.	0.019	09:57
05 Room Air	0.000	09:57
06 Std. Sol.	0.019	09:58
07 Room Air	0.000	09:58

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 777777

*DK*  
*Deb Kasher*

Operator Signature  
DEB KASHUR

Remarks:

*Low AC 0.021 AC*

Form 106-I8000

80-004189  
p. 114 of 120  
FM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = TOXL      8164.13.00 06/09  
04/01/2011      09:59

WET CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:59
02 Std. Sol.	0.294	10:00
03 Room Air	0.000	10:01
04 Std. Sol.	0.294	10:02
05 Room Air	0.000	10:02
06 Std. Sol.	0.295	10:03
07 Room Air	0.000	10:03

08 Sim Temp = 34.0°C

Simul Ser No = 0

Std Sol No = 0

County = 08

Oper No. = 777777

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

*Linearity test 0.294 AC*

Form 106-I8000

80-004189  
P. 115 of 120  
EM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004189  
Location = TOXL                8164.13.00 06/09  
04/01/2011                        09:45

Test	AC	Time
01 Room Air	0.000	09:46
02 *Subject Test	INT*	09:46
03 Room Air	0.000	09:47

\*Invalid Test  
Interferent Detected

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 05911080A1  
Cyl No = 17  
Expiration Date = 02/01/2013  
County = 08                          Oper No. = 777777

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

*Interferent test 0.10 AC + acetone*

Form 106-I8000

80-004189  
P. 116 of 120  
KTM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004189  
Location = TOXL                    8164.13.00 06/09  
04/01/2011                            09:51

Test	AC	Time
01 Diagnostic	OK	09:52
02 Room Air	RFI*	09:52
03 Room Air	0.000	09:52

\*Invalid Test  
Inhibited - RFI

Sub Name = DISCOVER, THE SPIRIT  
Sub DOB = 02/01/1992  
Sub Sex = Female                    Weight = 150  
Test = DUI                            Cit = NA  
Dr. Lic. = ND/DIS921456  
Lot No = 05911080A1  
Cyl No = 17  
Expiration Date = 02/01/2013  
County = 08                            Oper No. = 777777

I followed the Approved Method and the instructions displayed by the Intoxilyzer in conducting this test.

Deb Kasher

Operator Signature  
DEB KASHUR

Remarks:

RFI test

Form 106-I8000

80-004189  
P. 117 of 120  
DK

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer            Alcohol Analyzer  
North Dakota Model 8000        SN 80-004189  
Location = TOXL                8164.13.00 06/09  
04/01/2011                        09:30

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:30
02 Std. Gas	0.078	09:31
03 Room Air	0.000	09:31
04 Std. Gas	0.079	09:32
05 Room Air	0.000	09:32
06 Std. Gas	0.080	09:33
07 Room Air	0.000	09:33

Lot No = 05911080A1  
Cyl No = 17  
Exp Date = 02/01/2013  
County = 08                      Oper No. = 777777

Deb Kashur  
Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

80-004189  
P. 118 of 120  
JM

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = TOXL      8164.13.00 06/09  
04/01/2011      09:33

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:34
02 Std. Gas	0.079	09:34
03 Room Air	0.000	09:35
04 Std. Gas	0.080	09:35
05 Room Air	0.000	09:36
06 Std. Gas	0.079	09:36
07 Room Air	0.000	09:37

Lot No = 05911080A1  
Cyl No = 17  
Exp Date = 02/01/2013  
County = 08      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

80-004189  
P. 119 of 120  
AM



Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-004189  
Location = TOXL      8164.13.00 06/09  
04/01/2011      09:39

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	09:40
02 Std. Gas	0.079	09:40
03 Room Air	0.000	09:41
04 Std. Gas	0.080	09:41
05 Room Air	0.000	09:41
06 Std. Gas	0.079	09:42
07 Room Air	0.000	09:42

Lot No = 05911080A1  
Cyl No = 17  
Exp Date = 02/01/2013  
County = 08      Oper No. = 777777

Deb Kasher  
Operator Signature  
DEB KASHUR

Remarks:

*Calibration check 0.080 AC*

Form 106-I8000

80-004189  
p. 120 of 120  
AM