



**ETHANOL BREATH STANDARD CYLINDER REPORT**  
 OFFICE OF ATTORNEY GENERAL  
 CRIME LABORATORY DIVISION  
 Toxicology Section/Breath Alcohol Program  
 SFN 59282 (5-2011)

Chemical Test Operator Name (Print) <i>Shannon Liebold</i>		
Location <i>Mercer County Jail</i>	Intoxilyzer® Serial Number <i>80-003060</i>	
Gas Lot Number <i>1161108042</i>	Gas Cylinder Number <i>115</i>	Gas Expiration Date <i>06/01/2013</i>

Check When Done:

- 1. Scan/Enter Gas Cylinder Information (Level 1, Function S).
- 2. Perform an ACA Test (Level 1, Function C).  
 Write in the ACA Test Results (Reported to 3 Digits, Ex. 0.081).  
 A. 0. 079 AC  
 B. 0. 079 AC  
 C. 0. 079 AC
- 3. Complete the Top Portion of the Intoxilyzer® Record (SFN50496, Form 120-G) and Place it by the Intoxilyzer® for Use.
- 4. File Previous Intoxilyzer® Record (SFN50496, Form 120-G) at the Intoxilyzer® Location at the Agency.
- 5. Send the Following to the Crime Laboratory:
  - A. Ethanol Breath Standard Cylinder Report (SFN59282, Form 105-G).
  - B. ACA Test Record.

Chemical Test Operator Signature <i>Shannon Liebold</i>	Date <i>9-9-11</i>
Reviewed By (Crime Laboratory Use Only) <i>Deb Kashner</i>	Date <i>02 Apr 12</i>

*(Extra paperwork sent in)  
 PK*

Intoxilyzer Test Record and Checklist  
NDOAG Crime Lab. Div., Bismarck, ND 58501

CMI, Inc. Intoxilyzer      Alcohol Analyzer  
North Dakota Model 8000      SN 80-003060  
Location = STAT      8164.13.00 06/09  
09/09/2011      17:43

DRY CAL CHECK

Test	AC	Time
01 Room Air	0.000	17:44
02 Std. Gas	0.079	17:45
03 Room Air	0.000	17:45
04 Std. Gas	0.079	17:45
05 Room Air	0.000	17:46
06 Std. Gas	0.079	17:46
07 Room Air	0.000	17:47

Lot No = 11611080A2  
Cyl No = 115  
Exp Date = 06/01/2013  
County = 29      Oper No. = 020904



Operator Signature  
SHANNON P. LIEBOLD

Remarks: